

Electronic Signature Information

| | |
|----------------------|---|
| Name | DOC1620970 |
| Revision | 1 |
| Type | Controlled Document |
| Title | SN Quench Recovery Resistance Clarification |
| Originator | 212000253_michelle_g_sansbury |
| Release Date | |
| Obsolete Date | |

| Name | Reason For Change | File Size (Bytes) |
|---|---|-------------------|
| SN_DOC1620970_Magnet_Quench_Recovery_Resistance_R001.docx | SN_DOC1620970_Magnet_Quench_Recovery_Resistance_R001.docx | 363392 |

| Route | Signer | Function | Status | Comments | Completion Date |
|-----------|-------------------------------|------------------|---------|------------------------|--------------------------|
| R-6019691 | 212000253_michelle_g_sansbury | IB Test Engineer | Approve | Approved | 10 Nov 2014 13:32:42 GMT |
| R-6019691 | 212011747_timothy_j_havens | LSSD-LSI | Approve | Approved | 4 Nov 2014 20:54:43 GMT |
| R-6019691 | 212012451_juan_r_rosa | QARA | Review | Approve | 5 Nov 2014 20:48:24 GMT |
| R-6019691 | 212031211_yuri_ivovsky | LSD | Review | Reviewed and approved. | 7 Nov 2014 14:47:42 GMT |
| R-6019691 | 212062892_andrei_baikalov | Des and Dev | Approve | approve | 7 Nov 2014 19:47:30 GMT |

Periodic Review

There are no signatures or routes related to this business object.

Obsolescence Approval

There are no signatures or routes related to this business object.

- * Printed versions are For Reference Only *
- + Indicates a task was reassigned from an original assignee



Quench Recovery

Resistance Value Clarification



1. SERVICE NOTE CLASSIFICATION

- A GE Healthcare MR Class A eService Notes are copyright (c) protected to General Electric Company with all rights reserved. Duplication and distribution of any class of MR eService Notes in any form, by entities other than GE Healthcare, is strictly prohibited.
- C GE Healthcare MR Class C eService Notes are proprietary to General Electric Company. Disclosure to persons not employed by the General Electric Company or not in possession of an Advanced Service Agreement with GE Healthcare is prohibited. Duplication and distribution of Class C eService Notes in any form, by persons or entities other than by the General Electric Company, is strictly prohibited. Class C eService Notes are copyright (c) protected.
- M GE Healthcare MR Class M eService Notes are copyright (c) protected. GE Exclusive Use Only. Duplication and distribution of Class M eService Notes in any form, by persons or entities other than by the General Electric Company, is strictly prohibited.

2. PURPOSE

The purpose of this service note is to communicate a required update to the Quench Recovery section of 5495018, Magnet and Cryogen Manual for Passively Shimmed Magnets. Section 10.20.7 contains an incorrect reference to an open circuit.

3. SCOPE

This service note applies to

- A. 1.5T MR450w (HM Series) Magnets
- B. 3.0T MR750w (UA Series) Magnets

4. REFERENCES

- A. This service note is owned by Magnet and Gradient Engineering
- B. Questions pertaining to this Service Note should be directed to your Regional Magnet & Cryogen (MAC) Team Leader or the Magnet Support Team at the Online Center
- C. This service note will be incorporated into the next release of 5495018, Magnet and Cryogen Manual for Passively Shimmed Magnets.

5. SAFETY

Chapter 1 of the Magnet and Cryogens Subsystem Manual must be reviewed before the procedures in this Service Note are completed. All safety procedures must be adhered to at all times while performing this procedure.



Quench Recovery

Resistance Value Clarification



6. QUENCH RECOVERY SECTION PROCEDURE UPDATE
 Referencing section 10.20.7 (in 5495018, Revision 002) –
 PASSIVELY SHIMMED MAGNETS ONLY

- conformance with basic replacement.*
- 10.20.6. Purge the Recondenser with 4 bottles of 99.999% ("five nines") helium gas, inputting to the Recondenser shield vent line and exiting through the magnet's helium vent by opening Vent Valve V2. If applicable, use the DVw Deicing Kit (5342456) to assist with deicing the magnet. Refer to the Recondenser Deicing subsection of "Ice Removal Procedures."
 - 10.20.7. Connect a DVM between an engaged Main Lead Extension and an unpainted surface on the Cryostat to check for a Main Coil grounding fault. If an open circuit is not indicated, contact your Regional Magnet & Cryogen (MAC) Team Leader or the Magnet Support Team at the Online Center.
 - 10.20.8. Perform Magnet Electrical Checks as detailed within the Commissioning Check section of this manual.
 - 10.20.9. Check the Pressure Gauge for damage, replacing if required.
 - 10.20.10. Refill the magnet in conformance with "Liquid Helium Fill."
 - 10.20.11. Check Cryostat pressure after liquid helium fill and stabilization. Do not adjust Cryostat pressure at this time.

This section will be updated as follows:

- conformance with basic replacement.*
- 10.20.6. Purge the Recondenser with 4 bottles of 99.999% ("five nines") helium gas, inputting to the Recondenser shield vent line and exiting through the magnet's helium vent by opening Vent Valve V2. If applicable, use the DVw Deicing Kit (5342456) to assist with deicing the magnet. Refer to the Recondenser Deicing subsection of "Ice Removal Procedures."
 - 10.20.7. Connect a DVM between an engaged Main Lead Extension and an unpainted surface on the Cryostat to check for a Main Coil grounding fault. Expected resistance values are noted below; if readings outside these values are obtained, contact your Regional Magnet & Cryogen (MAC) Team Leader or the Magnet Support Team at the Online Center.
 - ◆ 1.5T MR450w (HM Series): 9,550 Ohms to 14,600 Ohms
 - ◆ 3.0T MR750w (UA Series): 6,600 Ohms to 9,000 Ohms
 - 10.20.8. Perform Magnet Electrical Checks as detailed within the Commissioning Check section of this manual.
 - 10.20.9. Check the Pressure Gauge for damage, replacing if required.
 - 10.20.10. Refill the maagnet in conformance with "Liquid Helium Fill."



Quench Recovery

Resistance Value Clarification



| Revision History | | | | | | |
|-------------------------------------|--------------------------|-------------------------------------|--|---|----------------------------|---------------------------------|
| REV | Affected Page(s) | Description of Content Changed | Reason for Change or Change Control Number | Author | Effective date of Document | Periodic review date |
| 01 | ALL | Initial Release | Error within 5495018. Copy error from legacy LCC manual. HM's and UA's will not have an open circuit - because the coil circuit grounding scheme on HM and UA is different from LCC which has open ground. | Michelle G. Sansbury | My Workshop Release | Does not Apply to Service Notes |
| Yes | No | NA | Enter Validation Information or Rationale if "NO" or "NA" for Validation or Training | | | |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Typo or Format Change | Copy error from legacy manual. | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Regulatory Filing | | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Validation | No change to existing process. Corrected resistance values. | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Training | None required. No change to existing process. | | |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Communication | Via email, upon MWS Release | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Labeling, (list) | | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Other Documents, (list) | | | |
| <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Other (describe) | | | |