


# TABLE OF CONTENTS

TABLE OF CONTENTS.....	1
1- Patient Information .....	2
2- Phantom Setting .....	5
3- Scan .....	6
4- Data Analysis .....	7
5- Data Sheet .....	11

Rev 2

### 1. Patient Information

1. Select Patient Information and Protocol Icon.



**PATIENT INFORMATION AND PROTOCOL ICON**  
ILLUSTRATION 1

2. Click [New Pt] button in "PATIENT REGISTER".



**[NEW PT] BUTTON**  
ILLUSTRATION 2

Rev 2

## 2. Scan Parameter Setting

1. Click **[New Pt]**.
2. Input the following data in "patient information".
  - Patient Data: geservice
  - Weight: 50(Kg)

*PATIENT INFORMATION*

Accession Number

Patient ID

Patient Name

Birth Date  Age  Sex

Weight  (Lb)  (Kg)

Rad  Refer

Req Number  Stastus

Description


History

**PATIENT INFORMATION**  
ILLUSTRATION 3

Rev 2

## 2. Scan Parameter Setting(continued)

3. Choose **GE**.
4. Click **Other**. Protocol window comes Up.
5. Select [**RF Coil SNR(A)**] from protocol.  
Select [**SNR Ax FOV25**] from series.
6. Verify that coil name shows **9-Inch GP**. If not, select from menu.
7. Select [**Accept**].



PATIENT INFORMATION  
ILLUSTRATION 1

8. Select [**Save Series**].





SAVE SERIES

Rev 2

### 3. Phantom Setting

1. Set the Sphere Phantom on the Phantom holder.
2. Set the Sphere Phantom and the Phantom holder in 9 inch GP Coil.
3. Set 2) on the table. See illustration about the phantom position.



**PHANTOM SETTING**  
ILLUSTRATION 4

4. Turn alignment light ON.
5. Advance the cradle to the position where the A-light beam hits at center of 9 inch GP Coil.
6. Press [**Landmark**] button.
7. Press [**Adv to Scan**] button to send the phantom into the Magnet center.

Rev 2

**4. Scan**

1. Select **Auto Prescan** button




**AUTO PRESCAN**  
ILLUSTRATION 5

2. Select **Manual Prescan** button.



**MANUAL PRESCAN**  
ILLUSTRATION 6

3. Click on "SCAN TR" and check that the projection is displayed.



**PROJECTION**  
ILLUSTRATION 7

Rev 2

### 4. Scan(continued)

- 4. Select **Scan** button



**SCAN**  
ILLUSTRATION 8

- 5. Verify that scanning has been started.

- 6. Select **Scan** button again.



**SCAN**  
ILLUSTRATION 9

- 7. Verify that scanning has been started.

### 5. Data Analysis


#### 5-1 Subtract Scan Images

- 1. Select Advantage Window Icon.



**AW ICON**  
ILLUSTRATION 10

- 2. Click **[Add/Sub]** button.




**ADD-SUB BUTTON**  
ILLUSTRATION 11

Rev 2

**5-1 Subtract Scan Images (continued)**

- 3. Subtract Image  
1st Image - 2nd Image

- 1. Click on "-".
- 2. Click on "**Accept Negative**
- 3. Select the Scan Image 1.
- 4. Click on "**Select Set**".(Left side)
- 5. Select the Scan Image 2.
- 6. Click "**Select Set**".(Right Side)
- 7. Click on "=".
- 8. Click "**Quit**".
- 9. The subtracted image is saved as **Series 100**\*\*



**HOW TO SUBTRACT**  
ILLUSTRATION 12

**5-2 Measure the Mean and SD of 1st image.**

- 1. Open the 1st Image by using the Miniviewer.
  - a. Select Advantage Window Icon.




**AW ICON**  
ILLUSTRATION 13

Rev 2


**5-2 Measure the Mean and SD of 1st image (continued)**

b. Select "Exam", "Series", and "Image" from the patient list and then press [MiniViewer].



**MINIVIEWER**  
ILLUSTRATION 14

2. From "Cursors and Measuring", select ellipse and set the area to (10000±500mm<sup>2</sup>)




**ELLIPTICAL CURSOR**  
ILLUSTRATION 15

3. Measure the Mean and SD.
4. Record the Data into the Data Sheet.
5. Press the keys at the same time to copy the current elliptical cursor :  
[Ctrl + Alt + F6] or [Ctrl + c].

Rev 2

**5-3 Measure the Mean and SD of 2nd image.**

1. Open the 2nd Image
  - a. Click on "Image +" button, and next image will come up.



**Next Image**  
ILLUSTRATION 16

2. Press the keys at the same time to paste the elliptical cursor:  
[Ctrl + Alt + F8] or [Ctrl + v].
3. Measure the Mean and SD.
4. Record the Data into the Data Sheet.

**5-4 Measure the Subtracted Image.**

1. Open the Subtracted Image.
  - a. Click on "Image +" button, and next image will come up.
2. Measure the Mean and SD of the Subtracted Image.
3. Record the Data into the Data Sheet.

**6. Data Sheet**

1. Record the Mean and SD values, and calculate the SNR:

$$SNR = \left( \frac{m1 + m2}{2} \times \sqrt{2} \right) / SD(s)$$

where, M1 is the 1st Image Mean Value and SD(s) is the Subtracted Image SD value

TABLE 6-1

**MEAN AND SD VALUES**

	MEAN	SD
1st Image		
2nd Image		
Subtracted Image	-----	

TABLE 6-2

**SNR VALUE**

Calculated SNR	SPEC
	≥ 189

Rev 2

**Revision History**

<b>Rev</b>	<b>Date</b>	<b>Author</b>	<b>Primary Reasons For Change</b>
0	Feb 9, 2001	Y.Masumo	Initial Releases
1	May 15, 2001	Y. Masumo	Coil setting and spec are updated.
2	Jan 28, 2002	Y. Masumo	Misc Correction.