

# TABLE OF CONTENTS

TABLE OF CONTENTS.....	1
1- Scan Parameter Setting .....	2
2- Phantom Setting .....	4
3- Scan .....	5
4- Data Analysis .....	7

Rev 2

### 1. Scan Parameter Setting

1. Click **[New Pt]**.
2. Input the following data in "patient information".
  - Patient Data: geservice
  - Weight: 50(Kg)

*PATIENT INFORMATION*

Accession Number

Patient ID

Patient Name

Birth Date  Age  Sex

Weight  (Lb)  (Kg)

Rad  Refer

Req Number  Stastus

Description


History

**PATIENT INFORMATION**  
ILLUSTRATION 1

Rev 2

1. Scan Parameter Setting (continued)

- 3. Choose GE.
- 4. Click **Other**. Protocol window comes Up.
- 5. Select [**System Check Head**] from protocol.  
Select [**Noise**] from series.
- 6. Select [**Accept**]. )



PATIENT INFORMATION  
ILLUSTRATION 2

- 7. Select Save Series.





SAVE SERIES  
ILLUSTRATION 3

Rev 2

## 2. Phantom Setting

1. Set cradle spacer to cradle.
2. Set head rest onto cradle.  
(Head Rest is used for positioning purpose only. Do Not set the phantom on the head coil.)
3. Set head coil to head rest..

**NOTE:**There are two type of table as following illustration.




**HEAD COIL SETTING 1**  
ILLUSTRATION 4

Rev 2

**2. Phantom Setting(Continued)**

- 4. Connect the coil connector box to table connector port.
- 5. Attach the coil cable to cable clamp of head rest.




**HEAD COIL SETTING 2**  
ILLUSTRATION 5

Rev 2

**2. Phantom Setting(Continued)**

6. Advance the cradle to the position where the A-light beam hits the center of coil.



**ALIGNMENT LIGHT**  
ILLUSTRATION 6

7. At front enclosure on scanner, press LANDMARK, then MOVE TO SCAN.

Rev 2


**3. Scan**

8. Select [Manual Prescan].



**MANUAL PRESCAN**  
ILLUSTRATION 7

9. Set TG=1, R1=13, R2=15. Then, set Center Frequency to “14.840000”.



**MANUAL PRESCAN**  
ILLUSTRATION 8

10. Select [Scan]..




**SCAN**  
ILLUSTRATION 9

11. Verify that scanning has been started.

Rev 2

### 3. Scan (continued)


12. Select [Manual Prescan].



#### MANUAL PRESCAN

ILLUSTRATION 10

13. Set TG=1, R1=13, R2=15. Then, set Center Frequency to "14.860000".



#### MANUAL PRESCAN

ILLUSTRATION 11

14. Select [Scan]..



#### SCAN

ILLUSTRATION 12

15. Verify that scanning has been started.

Rev 2


### 4. Data Analysis

1. Open the Image.

1. Select [Advantage Windows] Icon.




2. Select "Exam", "Series", and "Image" from the patient list and then press [Viewer].



**OPEN IMAGE**  
ILLUSTRATION 13

2. Measure Mean Value (Ms) of the brightest dot.

1. Set WW to 0.
2. Increase WL until only the one brightest dot can be seen on the Image.
3. Record the value of WL as Ms.





**DOT NOISE**  
ILLUSTRATION 14

Rev 2

**4. Data Analysis (continued)**

3. Measure Mean Value of the center of the Image(Mn) using circular cursor.

1. Select Elipse  from Graphic Menu.
2. Set the circular cursor to the center of the phantom. Be sure that there is no bright part inside of the circle.
3. ROI (Region of Interest) comes up to the right bottom of the Image automatically.
4. Record the SD and Mean value of ROI.



**BACKGROUND NOISE**  
ILLUSTRATION 15

4. Verify the following formula.

$$Ms \leq Mn \times 10$$

**Table 1:**

fc	Ms	Mn
14.840000		
14.860000		

Rev 2

**Revision History**

<b>Rev</b>	<b>Date</b>	<b>Auther</b>	<b>Primary Reasons For Change</b>
0	Jan 26, 2000	Y. Masumo	Initial Release
1	May 15, 2001	Y. Masumo	Head Coil Setting is updated.
2	Oct 22, 2001	Y. Masumo	Page 5: Added the purpose of Head Rest. Page 7 & Page 8: Corrected the center frequency values.