

# TABLE OF CONTENTS

<b>TABLE OF CONTENTS</b> .....	<b>1</b>
<b>1- Preparation</b> .....	<b>2</b>
1-1 Overview.....	2
1-2 Required Tool .....	2
1-3 Flowchart .....	3
1-4 Prerequisite .....	3
<b>2- Setup</b> .....	<b>4</b>
2-1 Mapping Fixture setup .....	4
2-2 Measurement Tool Setup .....	7
2-3 Probe Setup. ....	10
<b>3- Field Stabilization</b> .....	<b>11</b>
<b>4- Shim Data Gathering Program Setup</b> .....	<b>11</b>
4-1 Start Shim Data Gathering Program .....	11
4-2 Shim Data Gathering Program with Voice synthesizer .....	12
4-3 Shim Data Gathering Program without Voice synthesizer .....	18
4-4 Mechanical Plot Setting .....	22
<b>5- Data Gathering</b> .....	<b>24</b>
<b>6- Get the Result of Data Gathering</b> .....	<b>29</b>
6-1 Start the Program .....	29
6-2 Modify the Field Map File and Passive Shim .....	30
6-3 Mechanical or Passive Shimming .....	32
6-3-1 If $ Z2  \geq 15$ ppm... ..	33
6-3-2 If $ Z2  < 15$ ppm and $ Z1 ,  X ,  Y ,  ZX $ , or $ ZY  \geq 4$ ppm... ..	35
6-3-3 If $ Z2  < 15$ ppm and $ Z1 ,  X ,  Y ,  ZX $ , or $ ZY  < 4$ ppm... ..	39
<b>7- Passive Shimming</b> .....	<b>42</b>
7-1 Upper Shim Plate Shim Setting .....	43
7-2 Lower Shim Plate Shim Setting .....	46

Rev 4

**1. Preparation**

**1-1 Overview**

This procedure describes how to calibrate the magnet field distribution inside of the sphere.

There are two way for the shimming:

- a.Mechanical Shimming
- b.Passive Shimming

Mechanical Shimming adjust the  $Z^2$ , Z, X, Y, ZX, ZY using adjustment bolts attached to the Yoke of the Magnet. Passive shimming is performed if the pk-pk value is over 30ppm by attaching the shim plate on shim board. Shim plate position and number are calculated by software.

Passive shimming adjust the higher terms than  $Z^3$  and partial inhomogeneity


**1-2 Required Tool**

TABLE 1-1  
REQUIRED TOOL

Items	Check
T-wrench and shaft	
Tesla Meter (metro lab tesla meter PT2025 (or equivalent))	
Probe (No 3)	
Service Key	
Mapping Fixture	
Automated Shimming Data Collector (46-323012G1)(Option) (Voice Synthesizer (Option) and Remote Switch (Option) are included)	
Passive Shim Kit (N/ S checker, Shim chips)	

Rev 4

1-3 Flowchart



SHIMMING FLOWCHART  
ILLUSTRATION 1

1-4 Prerequisite


1. TCU Power is ON and magnet temperature is stabiliszed.
2. Drift Rate of Magnet Center freaquency must be below 50Hz/5minutes.  
(See "Section 3 Field Stabilization" for detail)

Rev 4

**2. Setup**

**2-1 Mapping Fixture setup**


1. Attach the two Beams to the Base Plate with flat headed screws.




**MAPPING FIXTURE SETUP 1**  
ILLUSTRATION 2

2. Attach Positioning Bar to the Base Plate.

Gear (or Shaft), Conn Plate, and Positioning Bar will be attached together. (Shaft is attached instead of gear at the upper side of the Positioning Bar)



The Guide Block, End Block, and Index Rod will come attached to the Base Plate (The gear is attached to the Positioning Bar). Attach O-Ring as shown.




**MAPPING FIXTURE SETUP 2**  
ILLUSTRATION 3

Rev 4

### 2-1 Mapping Fixture setup (continued)


3. Attach the Top Plate over the Beam and the Positioning Bar with flat headed screw.



**MAPPING FIXTURE SETUP 3**

ILLUSTRATION 4


4. Attach the two Base Plate under both sides of the Base Plate.



**MAPPING FIXTURE SETUP 5**

ILLUSTRATION 5

5. Install the Mapping Fixture with studs.





**MAPPING FIXTURE SETUP 6**

ILLUSTRATION 6

Rev 4

**2. Mapping Fixture setup (continued)**

- Adjust the gap between RF Shield and of the 4 corners of Mapping Fixture using spacers (1mm thickness or 2mm thickness).



Specification: GAP D - GAP A =  $\pm 1$ mm  
 GAP C - GAP B =  $\pm 1$ mm  
 GAP A - GAP B =  $\pm 1$ mm  
 GAP D - GAP C =  $\pm 1$ mm

**MAPPING FIXTURE SETUP 7**  
 ILLUSTRATION 7


Rev 4

## 2-1 Measurement Tool Setup

1. Turn system power OFF before cabling. Turn ON after cabling.
2. Set Up according to the following Illustration.


### Cable Setup (Ovev view) :

Refer to each detail setup followed after this illustration.



**CABLE SETUP (OVERVIEW)**  
ILLUSTRATION 8

Serial Port 8: To Tesla Meter  
Serial Port 7: To Voice Synthesizer (option)




**SERIAL CABLE**  
ILLUSTRATION 9

Rev 4

## 2-1 Measurement Tool Setup (continued)


### Tesla Meter Setup:

For METROLAB NMR Teslameter PT 2025



**TESLA METER SETUP**  
ILLUSTRATION 10

### Pre-Amp Setup:




**PREAMP BOX**  
ILLUSTRATION 11

Rev 4

## 2-1 Measurement Tool Setup (continued)

### Voice Synthesizer(Optional) Setup:


( OPTION)



**VOICE SYNTHESIZER (OPTION)**  
ILLUSTRATION 12

### Speaker(Optional) Setup:

( OPTION)




**SPEAKER (OPTION)**  
ILLUSTRATION 13

Rev 4


**2-2 Probe Setup.**

1. Place the Probe on the Probe Plate
2. Screw the Clamp over the Probe(4 screws).



**PROBE SETUP 1**  
ILLUSTRATION 14

3. Attach the Probe to the **Magnet physical center(Isocenter)** as shown. Match the number "1" of the Probe Plate to "7" of the Positioning Bar.



**PROBE SETUP 2**  
ILLUSTRATION 15

Rev 4


### 3. Field Stabilization

1. Verify that the probe is at the Magnet Isocenter.
2. Verify that the Physical Center Frequency is 14.85 MHz ± 10kHz.  
If not, adjust the temperature of the magnet and wait until the condition is satisfied.
3. Measure and record the frequency shown in Tesla Meter.
4. Wait 5 minutes.
5. Measure and record the frequency again.
6. If the difference between the 2 frequencies is:
  1. Below 50Hz/5minutes, go to Step 4 "Start Shim Data Gathering Program".
  2. Over 50Hz/5minutes, it is necessary to wait until the difference is below 50Hz/5minutes.

### 4. Shim Data Gathering Program Setup

#### 4-1 Start Shim Data Gathering Program

1. Set Service Key.
2. Select [**Service Desktop**] Icon.
3. Select [**Cal/Checks**] .
4. Select [**Shim Data Gathering**].
5. Select [**Start...**].



**SERVICE DESKTOP**  
ILLUSTRATION 16

6. Select your system configuration:
  - If you are using Voice synthesizer, Perform 4-2 “Shim Data Gathering Program with Voice synthesizer”.
  - If not, Perform 4-3 “Shim Data Gathering Program without Voice synthesizer”.

Rev 4

### 4-2 Shim Data Gathering Program with Voice synthesizer

Output/Prompts	Input/Comments
<pre> **** MR Magnet Mechanical Shim Data Collection **** Do you use the vocie synthesizer ? [1:Yes / 2:No] --&gt; 1 MR Magnet Mechanical Shim Data Collection Enter number of selection: 1 - Normal operation 2 - Run mechanical plot only 3 - Calculate shim currents from existing data 4 - Display plot data 5 - Display current data 6 - Display magnet file information 7 - Misc functions 8 - Drift Test 9 - Escape to shell 10 - Exit --&gt; 7 Miscellaneous Functions Enter number of selection: 1 - System checkout 2 - Synthesizer volume adjust 3 - Teslameter bandwidth 4 - Teslameter # samples to average/point 5 - Teslameter retry count 6 - Select/Deselect Switch 7 - Enable/Disable switch feedback 8 - Display plot format (volume map) file information 9 - Display current plot format data 10 - Exit --&gt; 1 System Checkout Enter number of selection: 1. Checkout voice synthesizer 2. Checkout teslameter 3. Checkout remote switch operation 4. Exit --&gt; 1 ##### -----                     </pre>	<pre> 1 [Return]  7 [Return]  1 [Return]  1 [Return]                     </pre>



Rev 4

4-2 Shim Data Gathering Program with Voice synthesizer (continued)

Output/Prompts	Input/Comments
<pre> Enter number of selection:   1. Checkout voice synthesizer   2. Checkout teslameter   3. Checkout remote switch operation   4. Exit --&gt; 2       METROLAB TESLAMETER CHECKOUT ----- Press space bar to attempt communication, return to exit <b>(If connection is correct, the following message is read:)</b>   Teslameter is communicating <b>(If not:)</b>   No data received - check cable and connection   Press space bar to attempt communication, return to exit   System Checkout   Enter number of selection:     1. Checkout voice synthesizer     2. Checkout teslameter     3. Checkout remote switch operation     4. Exit --&gt; 4       Miscellaneous Functions   Enter number of selection:     1 - System checkout     2 - Synthesizer volume adjust     3 - Teslameter bandwidth     4 - Teslameter # samples to average/point     5 - Teslameter retry count     6 - Select/Deselect Switch     7 - Enable/Disable switch feedback     8 - Display plot format (volume map) file information     9 - Display current plot format data     10 - Exit --&gt; 2                     </pre>	<pre> 4 [Return]  Hit Space bar  Check connection and hit space bar [Return]  4 [Return]  2 [Return]                     </pre>



Rev 4

4-2 Shim Data Gathering Program with Voice synthesizer (continued)

Output/Prompts	Input/Comments
<pre> Miscellaneous Functions   Enter number of selection:     1 - System checkout     2 - Synthesizer volume adjust     3 - Teslameter bandwidth     4 - Teslameter # samples to average/point     5 - Teslameter retry count     6 - Select/Deselect Switch     7 - Enable/Disable switch feedback     8 - Display plot format (volume map) file information     9 - Display current plot format data    10 - Exit   --&gt; 5       Teslameter Retry Count Adjustment This will allow you to change the number of retries... Current retry count is 10 Enter new retry count (2-10):10 Miscellaneous Functions   Enter number of selection:     1 - System checkout     2 - Synthesizer volume adjust     3 - Teslameter bandwidth     4 - Teslameter # samples to average/point     5 - Teslameter retry count     6 - Select/Deselect Switch     7 - Enable/Disable switch feedback     8 - Display plot format (volume map) file information     9 - Display current plot format data    10 - Exit   --&gt; 6       Select/Deselect Switch You can change your decision about whether or not to use the remote switch. ... Select/Deselect Switch:Y                     </pre>	<pre> 5 [Return] 10 [Return] 6 [Return] Y [Return]                     </pre>

Rev 4

**4-2 Shim Data Gathering Program with Voice synthesizer (continued)**

Output/Prompts	Input/Comments
<pre> Miscellaneous Functions Enter number of selection:   1 - System checkout   2 - Synthesizer volume adjust   3 - Teslameter bandwidth   4 - Teslameter # samples to average/point   5 - Teslameter retry count   6 - Select/Deselect Switch   7 - Enable/Disable switch feedback   8 - Display plot format (volume map) file information   9 - Display current plot format data  10 - Exit --&gt; 10                     </pre>	<p>10 [Return]</p>

Rev 4

4-3 Shim Data Gathering Program without Voice synthesizer

Output/Prompts	Input/Comments
<pre> Do you use the voice synthesizer ? [1:Yes / 2:No] --&gt; 2 MR Magnet Mechanical Shim Data Collection Enter number of selection: 1 - Normal operation 2 - Run mechanical plot only 3 - Calculate shim currents from existing data 4 - Display plot data 5 - Display current data 6 - Display magnet file information 7 - Misc functions 8 - Drift Test 9 - Escape to shell 10 - Exit --&gt; 7 Miscellaneous Functions Enter number of selection: 1 - System checkout 2 - Synthesizer volume adjust 3 - Teslameter bandwidth 4 - Teslameter # samples to average/point 5 - Teslameter retry count 6 - Select/Deselect Switch 7 - Enable/Disable switch feedback 8 - Display plot format (volume map) file information 9 - Display current plot format data 10 - Exit --&gt; 1 System Checkout Enter number of selection: 1. Checkout voice synthesizer 2. Checkout teslameter 3. Checkout remote switch operation 4. Exit --&gt; 2                     </pre>	<pre> 2 [Return]  7 [Return]  1 [Return]  2 [Return]                     </pre>

Rev 4

4-3 Shim Data Gathering Program without Voice synthesizer (continued)

Output/Prompts	Input/Comments
<pre> **** MR Magnet Mechanical Shim Data Collection ****           METROLAB TESLAMETER CHECKOUT           -----           Press space bar to attempt communication, return to exit           (If connection is correct, the following message is read:)           Teslameter is communicating           (If not:)           No data received - check cable and connection            Press space bar to attempt communication, return to exit           System Checkout           Enter number of selection:           1. Checkout voice synthesizer           2. Checkout teslameter           3. Checkout remote switch operation           4. Exit           --&gt; 4            Miscellaneous Functions           Enter number of selection:           1 - System checkout           2 - Synthesizer volume adjust           3 - Teslameter bandwidth           4 - Teslameter # samples to average/point           5 - Teslameter retry count           6 - Select/Deselect Switch           7 - Enable/Disable switch feedback           8 - Display plot format (volume map) file information           9 - Display current plot format data           10 - Exit           --&gt; 3           </pre>	<pre> Hit Space bar           (If connection is correct)            Check connetion and hit           space bar           [Return]            4 [Return]            3 [Return]           </pre>

Rev 4

4-3 Shim Data Gathering Program without Voice synthesizer (continued)

Output/Prompts	Input/Comments
<p>Teslameter Bandwidth Adjustment</p> <p>This will allow you to adjust the bandwidth, or in other words the maximum frequency deviation allowed among the samples for a given point before the data is accepted. Once accepted, the samples are averaged for a given point.</p> <p>...</p> <p>Current bandwidth is 16 hz</p> <p>Enter new bandwidth (3-30 Hz):30</p> <p>Miscellaneous Functions</p> <p>Enter number of selection:</p> <ul style="list-style-type: none"> <li>1 - System checkout</li> <li>2 - Synthesizer volume adjust</li> <li>3 - Teslameter bandwidth</li> <li>4 - Teslameter # samples to average/point</li> <li>5 - Teslameter retry count</li> <li>6 - Select/Deselect Switch</li> <li>7 - Enable/Disable switch feedback</li> <li>8 - Display plot format (volume map) file information</li> <li>9 - Display current plot format data</li> <li>10 - Exit</li> </ul> <p>--&gt; 4</p> <p>Number of Readings to Average per Point</p> <p>This will allow you to adjust the number of readings that will be averaged for each point. The default...</p> <p>Current sample count is 3 readings</p> <p>Enter new sample count (2-10):3</p> <p>Miscellaneous Functions</p> <p>Enter number of selection:</p> <ul style="list-style-type: none"> <li>1 - System checkout</li> <li>2 - Synthesizer volume adjust</li> <li>3 - Teslameter bandwidth</li> <li>4 - Teslameter # samples to average/point</li> <li>5 - Teslameter retry count</li> <li>6 - Select/Deselect Switch</li> <li>7 - Enable/Disable switch feedback</li> <li>8 - Display plot format (volume map) file information</li> <li>9 - Display current plot format data</li> <li>10 - Exit</li> </ul> <p>--&gt; 5</p>	<p>30 [Return]</p> <p>4 [Return]</p> <p>3 [Return]</p> <p>5 [Return]</p>

Rev 4

4-3 Shim Data Gathering Program without Voice synthesizer (continued)

Output/Prompts	Input/Comments
<pre> This will allow you to change the number of retries... Current retry count is 10 Enter new retry count (2-10):10 Miscellaneous Functions Enter number of selection: 1 - System checkout 2 - Synthesizer volume adjust 3 - Teslameter bandwidth 4 - Teslameter # samples to average/point 5 - Teslameter retry count 6 - Select/Deselect Switch 7 - Enable/Disable switch feedback 8 - Display plot format (volume map) file information 9 - Display current plot format data 10 - Exit --&gt; 6 Select/Deselect Switch You can change your decision about whether or not to use the remote switch. ... Select/Deselect Switch:Y Miscellaneous Functions Enter number of selection: 1 - System checkout 2 - Synthesizer volume adjust 3 - Teslameter bandwidth 4 - Teslameter # samples to average/point 5 - Teslameter retry count 6 - Select/Deselect Switch 7 - Enable/Disable switch feedback 8 - Display plot format (volume map) file information 9 - Display current plot format data 10 - Exit --&gt; 10                     </pre>	<pre> 10 [Return]  6 [Return]  Y [Return]  10 [Return]                     </pre>

Rev 4

### 4-4 Mechanical Plot Setting

1. Perform as follows:

Output/Prompts	Input/Comments
<pre>MR Magnet Mechanical Shim Data Collection Enter number of selection: 1 - Normal operation 2 - Run mechanical plot only 3 - Calculate shim currents from existing data 4 - Display plot data 5 - Display current data 6 - Display magnet file information 7 - Misc functions 8 - Drift Test 9 - Escape to shell 10 - Exit --&gt; 2  ---(Message comes here)--- Press Enter threshold value:</pre>	<pre>2 [Return]  [Return]</pre>

2. After selecting #2, enter the following data.:

Output/Prompts	Input/Comments
<pre>Enter magnet serial number: YC*** [Return] Is this GE or Oxford magnet ( enter GE, Oxford, or press RETURN to Default to GE: [Return] Enter magnet location;(1 to 15 chars): ***** [Return] Enter magnet type (or press, RETRUN for type s2): ASM Is the magnet ramped normal (red to red, black to black) (y/n):Y Enter volume type (c6,s40)(or press RETURN for s40):pmh [return]</pre>	<pre>Input serial number [Return] Input Hospital Name ASM [Return] Y pmh [Return]</pre>

Rev 4

**4-4 Mechanical Plot Setting (Continued)**

3. Hit **[Return]** key or **<Remote Switch>** when the information is correct. Enter "i", "d", or "s" if needed. (Refer to the following.):

Output/Prompts	Input/Comments
<pre>MR Magnet Mechanical Plotting   This iteration will be for plot a   GE Magnet   Magnet serial number: YC***   Magnet location      : *****   Magnet type : ASM   Volume type : pmh   Magnet IS ramped normal   Press &lt;return&gt; to continue   Enter i to re-run previous iteration   Enter d to alter magnet data for serial # **   Enter s to select a new serial number   --&gt;<b>[Return]</b></pre>	<p><b>[Return]</b></p>


**PLOTTING DATA STARTS RUNNNING.**

Output/Prompts	Input/Comments
<pre>GE Magnet Mechanical Plotting   This iteration will be for plot v   Shim data collection in progress, press any key to pause...   Current probe position: Isocenter</pre>	

Rev 4

**5. Data Gathering**

1. Confirm that the probe is at magnet isocenter.




**PROBE SETUP AT MAGNET CENTER**  
ILLUSTRATION 17


If you are using Voice synthesizer, press remote switch.

If not, press **[Return]** key on keyboard. Plotting Data starts running.

2. Set the probe to "Plane 1"



1. Remove Stopper and Pin.
2. Move the Probe to "1".
3. Fix the Probe to Positioning Bar with Stopper and Pin.



**PROBE SETUP AT PLANE 1**  
ILLUSTRATION 18



If you are using Voice synthesizer, press remote switch.

If not, press **[Return]** key on keyboard. Plotting Data starts running.

Rev 4


**5. Data Gathering (continued)**

3. Set the probe to next plane.




**For Plane1 to Plane7**

1. Remove Stopper and Pin.
2. Move the Probe
3. Fix the Probe to Positioning Bar with Stopper and Pin.



**For Plane8 to Plane12**

1. Remove Stopper and Pin.
2. Move the Probe.
3. Fix the Probe to Positioning Bar with Stopper and Pin.




**NEXT PLANE**  
ILLUSTRATION 19

Rev 4

**5. Data Gathering (continued)**

4. Measure 0 degree thru 330 degree, step 30 degree..



**ROTATION**  
ILLUSTRATION 20


If you are using Voice synthesizer, press remote switch.  
If not, press **[Return]** key on keyboard. Plotting Data starts running.

5. Go back to step 3 until Plane12 is performed.


Rev 4

**5. Data Gathering (continued)**

6. Set the probe to "Plane 13"



1. Remove Stopper and Pin.
2. Move the Probe to 13.
3. Fix the Probe to Positioning Bar with Stopper and Pin.




**PROBE SETUP AT PLANE 13**  
ILLUSTRATION 21

If you are using Voice synthesizer, press remote switch.  
If not, press **[Return]** key on keyboard. Plotting Data starts running.

Rev 4

5. Data Gathering (continued)

7. Measure the frequency at Magnet Isocenter.



PROBE SETTING AT MAGNET CENTER  
ILLUSTRATION 22


If you are using Voice synthesizer, press remote switch.  
If not, press **[Return]** key on keyboard. Plotting Data starts running.

Rev 4

## 6. Get the Result of Data Gathering

### 6-1 Start the Program

1. Set Service Key.
2. Select [**Service Desktop**] Icon.
3. Select [**Cal/Checks**] .
4. Select [**Passive Shimming**].
5. Select [**Start...**].



**SERVICE DESKTOP**  
ILLUSTRATION 23

6. The following menu will be displayed.

Output/Prompts	Input/Comments
<pre> &lt;&lt;&lt; Permanent Magnet Shimming MENU &gt;&gt;&gt;  1. Field Map File :           gemech_SHIM-C_v.plot 2. First Passive? :          yes 0. Accept Defaults and Run Analysis. q. Exit Passive/Mech Shim  Enter the index number to change the default: (0..2) [0] :                     </pre>	



Rev 4

### 6-2 Modify the Field Map File and Passive Shim

3. After the Last Passive Shim File was selected, Menu is displayed as follows.

Output/Prompts	Input/ Comments
<pre> &lt;&lt;&lt; Permanent Magnet Shimming MENU &gt;&gt;&gt;    1. Field Map File :  gemech_YC008_v.plot   2. First Passive? :  no       Last Passive Shim File :  <u>shim_loc.gemech_YC008_a.plot</u>    0. Accept Defaults and Run Analysis.   q. Exit Passive/Mech Shim  Enter the index number to change the default:  (0..2) [0] :</pre>	<pre> &lt;- Selected Shim File</pre>

Rev 4

### 6-3 Mechanical or Passive Shimming

4. After all setting is done, enter "0" and hit [Return] key.

Output/Prompts	Input/Comments
<pre> &lt;&lt;&lt; Permanent Magnet Shimming MENU &gt;&gt;&gt;    1. Field Map File :                gemech_SHIM-C_v.plot   2. First Passive? :                no   0. Accept Defaults and Run Analysis.   q. Exit Passive/Mech Shim  Enter the index number to change the default: (0..12) [0] :                     </pre>	<p><b>0 [Return]</b></p>

5. Result of gathered data will be displayed.

Rev 4

6-3-1 If  $|Z2| \geq 15$  ppm...

1. If  $|Z2| \geq 15$  ppm, following is displayed.

Output/Prompts	Input/Comments
<pre> Writing to setup file Writing to init file Executing pmshim_ssm  Isocenter Frequency      : Start      = 14844273 [Hz]                           : End          = 14844276 [Hz]                           : Difference =      3 [Hz]  (Field Map File): gemech_SHIM-C_v.plot  &lt;&lt;&lt; Display Current Field Map File [difference in unit of PPM] &gt;&gt;&gt;  deg  ( 1) ( 2) ( 3) ( 4) ( 5) ( 6) ( 7) ( 8) ( 9) (10) (11) (12) (13) --- -----   0  -25. -24. -11.  0. -4.  11.  14.  20.  3.  3. -2. -11. -9.  30  -25. -15. -11. -9.  4.  14.  9.  17.  6. -11. -1. -7. -9.     : 330  -25. -21. -11. -10.  3.  11.  21.  10.  3. -2. -1. -6. -9. --- ----- avg  -25. -17. -15. -6.  2.  14.  17.  13.  5. -1. -4. -6. -9. MxMn   0.  14.  14.  13.  17.  8.  23.  15.  14.  23.  14.  9.  0.  Peak to Peak = 53.9 PPM  Axial coefficients ...   z 1=  -4.3   <u>z 2= -17.6</u>  Transverse coefficients ...   x=  -0.4   y=  -0.7   zx=  0.1   zy= -0.2  ----- &lt;&lt;&lt; Mech_Shim Output Summary &gt;&gt;&gt; -----  Z2 Mech_shim is recommended.  <u>Shift center freq to 615 [Hz] by turning upper handle</u> <u>Shift center freq to 615 [Hz] by turning lower handle</u>  Write Shim Files ? (N,Y) [Y]: </pre>	<pre> Start:CF Value at the start time of Gathering End: CF Value at the last time of Gathering Difference: Start - End. It should be 450Hz or lower.  If the difference is 450 or higer. Recheck that the Drift rate is within 50Hz/5min. Perform the gathering again within 30~40 min.  If  Z2  value is higher than 15ppm, it is neces- sary to adjust Z2.  Target Frequency = CF + the value here  Y [Return] to save the output file N [Return] If you calcu- late as trial, choose "N". Or, the result of this calcuration influence to the next calculation. </pre>

Rev 4

**6-3-1f  $|Z2| \geq 15$  ppm... (continued)****Note**


To avoid large center frequency drift, you have to perform mechanical shimming as soon as you can after calculation.

**Note**

Remove the cylindrical shaped insulator covering the mechanical shim point when performing the mechanical shimming. Do not remove another insulator. After the mechanical shimming is finished and restoring cylindrical shaped insulator, be careful not to push the insulator too much. Otherwise, it will be hard to remove it

1. Turn the upper adjustment bolt until Target Frequency is displayed on Tesla Meter.
  - a. Read Center Frequency of displayed on Tesla Meter.
  - b. Calculate the Target Frequency

Target Frequency = Center Frequency + Value of Mech Shim Output Summary




Turn the upper adjustment bolt using T-wrench

**Z2 ADJUSTMENT (UPPER)**

ILLUSTRATION 24

2. Turn the lower adjustment bolt until Target Frequency is displayed on Tesla Meter.
  - a. Read Center Frequency of displayed on Tesla Meter.  
(The value after the upper adjustment)
  - b. Calculate the Target Frequency

Target Frequency = Center Frequency + Value of Mech Shim Output Summary



Turn the lower adjustment bolt using T-wrench

**Z2 ADJUSTMENT (LOWER)**

ILLUSTRATION 25

Rev 4

6-3-2 If  $|Z2| < 15$  ppm and  $|Z1|, |X|, |Y|, |ZX|,$  or  $|ZY| \geq 4$  ppm...

1. If  $|Z2| < 15$  ppm and  $|Z1|, |X|, |Y|, |ZX|,$  or  $|ZY| \geq 4$  ppm, following is displayed.

Output/Prompts	Input/Comments
Writing to setup file Writing to init file Executing pmshim_ssm	
Isocenter Frequency : Start = 14844273 [Hz] : End = 14844276 [Hz] : Difference = 3 [Hz]	<b>Start:</b> CF Value at the atart time of Gathering <b>End:</b> CF Value at the last time of Gathering <b>Difference:</b> Start - End.
(Field Map File): gemech_SHIM-C_v.plot	It should be 450Hz or lower.
<<< Display Current Field Map File [difference in unit of PPM] >>>	If the difference is 450 or higer. Recheck that
<pre> deg  ( 1) ( 2) ( 3) ( 4) ( 5) ( 6) ( 7) ( 8) ( 9) (10) (11) (12) (13) --- -----   0  -25. -24. -11.  0. -4.  11.  14.  20.  3.  3. -2. -11. -9.  30  -25. -15. -11. -9.  4.  14.  9.  17.  6. -11. -1. -7. -9.       330  -25. -21. -11. -10.  3.  11.  21.  10.  3. -2. -1. -6. -9. --- ----- avg  -25. -17. -15. -6.  2.  14.  17.  13.  5. -1. -4. -6. -9. MxMn   0.  14.  14.  13.  17.  8.  23.  15.  14.  23.  14.  9.  0.                     </pre>	the Drift rate is within 50Hz/5min. Perform the gathering again within 30~40 min.
Peak to Peak = 53.9 PPM Axial coefficients ... z 1= -4.3 z 2= -14.6 Transverse coefficients ... x= -0.4 y= -0.7 zx= 0.1 zy= -0.2	If $ Z1 $ value is higher than 4ppm, it is neces- sary to adjust Z1 by mechanical shimming.
<<< Mech_Shim Output Summary >>>	
Z1 Mech_shim is recommended.	
Shift center freq to 204 [Hz] by turning Top handle	
Shift center freq to -204 [Hz] by turning Bottom handle	
Trans Mech_shim is recommended.	<b>Target Frequency</b> = CF + the value here
TA= 2.0 (CW ) TB= 1.5 (CCW) TC= 2.0 (CCW) TD= 1.5 (CW )	Value of rotation to adjust the Transverse Term.
BA= 2.0 (CW ) BB= 1.5 (CW ) BC= 2.0 (CCW) BD= 1.5 (CCW)	Y [Return] to save the output file
Write Shim Files ? (N,Y) [Y]:	N [Return] If you calcu- late as trial, choose "N". Or, the result of this calculation influence to the next calculation.

Rev 4

**6-3-2If  $|Z2| < 15$  ppm and  $|Z1|, |X|, |Y|, |ZX|,$  or  $|ZY| \geq 4$  ppm... (continued)****Note**

To avoid large center frequency drift, you have to perform mechanical shimming as soon as you can after calculation.

1. Turn the upper adjustment bolt until Target Frequency is displayed on Tesla Meter.
  - a. Read Center Frequency of displayed on Tesla Meter.
  - b. Calculate the Target Frequency


Target Frequency = Center Frequency + Value of Mech Shim Output Summary

**Note**

Remove the cylindrical shaped insulator covering the mechanical shim point when performing the mechanical shimming. Do not remove another insulator.

**CAUTION**

**For Z1 Adjustment, there is a mechanical stopper at adjustment screw. If you feel screw is reached to the stopper, DO NOT turn the screw.**




Turn the upper adjustment bolt using T-wrench

**Z1 ADJUSTMENT (UPPER)**

ILLUSTRATION 26

2. Turn the lower adjustment bolt until Target Frequency is displayed on Tesla Meter.
  - a. Read Center Frequency of displayed on Tesla Meter.  
(The value after the upper adjustment)
  - b. Calculate the Target Frequency

Target Frequency = Center Frequency + Value of Mech Shim Output Summary



Turn the lower adjustment bolt using T-wrench

**Z1 ADJUSTMENT (LOWER)**

ILLUSTRATION 27

Rev 4


**6-3-2f  $|Z2| < 15$  ppm and  $|Z1|, |X|, |Y|, |ZX|$ , or  $|ZY| \geq 4$  ppm... (continued)**

3. Turn the adjustment bolts at TA, TB, TC, TD, BA, BB, BC, BD according to the result file.

**Note**


Remove the cylindrical shaped insulator covering the mechanical shim point when performing the mechanical shimming. Do not remove another insulator.

After the mechanical shimming is finished and restoring cylindrical shaped insulator, be careful not to push the insulator too much. Otherwise, it will be hard to remove it.



**Note**

Be careful when turning the adjusting bolts of Top Yoke. Look up and turn the bolt to the direction described in the result file. See following illustration.



CW and CCW Direction Of Top Yoke

**MECHANICAL SHIMMING  
ILLUSTRATION 28**

Rev 4

6-3-2If  $|Z2| < 15$  ppm and  $|Z1|, |X|, |Y|, |ZX|, \text{ or } |ZY| \geq 4$  ppm... (continued)



If following message is appeared at the monitor, do not turn the screw.  
Mechanical Shimming cannot be performed any more.  
Perform passive shimming according to the instruction displayed at the monitor.

**\*\*CAUTION! Screws must not be turned any more.\*\***

Hit RETURN Key, and then start the passive-shim.

Rev 4

6-3-3 If |Z2| < 15 ppm and |Z1|, |X|, |Y|, |ZX|, or |ZY| < 4 ppm...

1. If |Z2| < 15 ppm and |Z1|, |X|, |Y|, |ZX|, or |ZY| < 4 ppm, following is displayed.

Output/Prompts	Input/Comments
<pre> Writing to setup file Writing to init file Executing pmshim_ssm  Isocenter Frequency      : Start      = 14844273 [Hz]                           : End        = 14844276 [Hz]                           : Difference =      3 [Hz] (Field Map File): gemech_SHIM-C_v.plot  &lt;&lt;&lt; Display Current Field Map File [difference in unit of PPM] &gt;&gt;&gt; deg  ( 1) ( 2) ( 3) ( 4) ( 5) ( 6) ( 7) ( 8) ( 9) (10) (11) (12) (13) --- -----   0  -25. -24. -11.  0. -4.  11.  14.  20.  3.  3. -2. -11. -9.  30  -25. -15. -11. -9.  4.  14.  9.  17.  6. -11. -1. -7. -9.       : 330  -25. -21. -11. -10.  3.  11.  21.  10.  3. -2. -1. -6. -9. --- ----- avg  -25. -17. -15. -6.  2.  14.  17.  13.  5. -1. -4. -6. -9. MxMn   0.  14.  14.  13.  17.  8.  23.  15.  14.  23.  14.  9.  0.  Peak to Peak = 53.9 PPM Axial coefficients ...   z 1=  -4.3   z 2= -14.6 Transverse coefficients ...   x=  -0.4   y=  -0.7   zx=  0.1   zy= -0.2  ----- &lt;&lt;&lt; Mech_Shim Output Summary &gt;&gt;&gt; ----- No need mechanical shimmming.  ----- &lt;&lt;&lt; Passive Shim Output Summary &gt;&gt;&gt; ----- Measured homogeneity in PPM  : 65.8 Calculated PPM (LP Unrounded) : 16.9 Calculated PPM (LP Rounded)  : 51.9 Calculated PPM (DS added)    : 32.4       : Location    U1  U2  U3  U4  U5  U6  U7   L1  L2  L3  L4  L5  L6  L7 ----- ----- Shim (mm)   4   4   7   7   7  12  12   4   4   7   7   7  12  12   SUM ----- ----- Current     0   0   0   0   0   0   0   0   0   0   0   0   0   0    0 Delta       0   3   4   4   1   4   1   0   0   0   0   0   1  11   29 ----- ----- Total       0   3   4   4   1   4   1   0   0   0   0   0   1  11   29 ----- ----- </pre>	<pre> &lt;- If measured homogeneity is under 30ppm and it is not necessary to perform passive shimmming, record this value in data sheet as a final center frequency.  If mechanical shim is not needed, this message is displayed.  If Measured Homogeneity is over 30ppm, Passive Shim Program will run.  LP:Linear Program Method DS:Direct Serch Method  &lt;- Current shim mag No. &lt;- Shim Mag No.to be added </pre>

Rev 4

6-3-3If |Z2| < 15 ppm and |Z1|, |X|, |Y|, |ZX|, or |ZY| < 4 ppm... (continued)

Output/Prompts													Input/Comments
Display Shim Direction and Shim Location *** TOP DELTA ***													Delta
INDEX	Rotation (Degrees)												
(No.)	0	30	60	90	120	150	180	210	240	270	300	330	
U1	0	0	0	0	0	0	0	0	0	0	0	0	
U2	0	0	N1	0	0	0	0	N1	0	0	S1	0	
U3	0	0	0	0	0	N1	0	0	0	N1	N1	N1	
U4	0	0	0	0	0	S1	0	0	0	S1	S1	S1	
U5	0	0	0	0	0	0	0	0	0	0	0	N1	
U6	0	0	0	0	N1	N1	-1	0	0	N1	0	0	
U7	0	0	-2N1	0	0	0	0	0	0	0	S1	0	
Display Shim Direction and Shim Location *** BOTTOM DELTA ***													See Next Page how to read this table.
INDEX	Rotation (Degrees)												
(No.)	0	30	60	90	120	150	180	210	240	270	300	330	
L1	0	0	0	0	0	0	0	0	0	0	0	0	
L2	0	0	0	0	0	0	0	0	0	0	0	0	
L3	0	0	0	0	0	0	0	0	0	0	0	0	
L4	0	0	0	0	0	0	0	0	0	0	0	0	
L5	0	0	0	0	0	0	0	0	0	0	0	0	
L6	0	0	S1	0	0	0	0	0	0	0	0	0	
L7	N2	0	0	-3N1	0	0	N2	N1	N3	0	0	0	
Display Shim Direction and Shim Location *** TOP TOTAL ***													Total
INDEX	Rotation (Degrees)												
(No.)	0	30	60	90	120	150	180	210	240	270	300	330	
U1	0	0	0	0	0	0	0	0	0	0	0	0	
U2	0	0	N1	0	0	0	0	N1	0	0	S1	0	
U3	0	0	0	0	0	N1	0	0	0	N1	N1	N1	
U4	S1	0	0	0	0	S1	0	0	0	S2	S1	S1	
U5	0	0	0	0	N1	0	0	0	0	0	0	N1	
U6	0	0	0	0	N1	N1	0	0	0	N1	0	0	
U7	0	0	0	0	0	0	0	0	0	0	S1	0	
Display Shim Direction and Shim Location *** BOTTOM TOTAL ***													Y [Return] to save the output file N [Return] not to save the output file
INDEX	Rotation (Degrees)												
(No.)	0	30	60	90	120	150	180	210	240	270	300	330	
L1	0	0	0	0	0	0	0	0	0	0	0	0	
L2	0	0	0	0	0	0	0	0	0	0	0	0	
L3	0	0	0	0	0	0	0	0	0	0	0	0	
L4	0	0	0	0	0	0	0	N1	0	N1	0	0	
L5	0	0	0	N1	0	0	N1	0	0	0	N1	N1	
L6	0	0	S2	0	0	0	0	0	0	0	0	0	
L7	N2	0	0	N1	S1	0	N2	N1	N3	0	0	0	
Write Shim Files ? (Y,N) [Y]:													
y													
-----Output file-----													
shim_loc.gemech_YC3_e.plot													
-----													
Press [Enter] to quit -->													

Rev 4

6-3-3If |Z2| < 15 ppm and |Z1|, |X|, |Y|, |ZX|, or |ZY| < 4 ppm... (continued)

**Note**

Refer to the following illustration to read the table.

Remove two Shim Magnets and add one "N" Shim Magnet

Remove one Shim Magnet

Add one "N" Shim Magnet

**For Top Shim Plate**

INDEX (No.)	Display Shim Direction and Shim Location											
	0	30	60	90	120	150	180	210	240	270	300	330
U1	0	0	0	0	0	0	0	0	0	0	0	0
U2	0	0	N1	0	0	0	0	N1	0	0	S1	0
U3	0	0	0	0	0	N1	0	0	0	N1	N1	N1
U4	0	0	0	0	0	S1	0	0	0	S1	S1	S1
U5	0	0	0	0	0	0	0	0	0	0	0	N1
U6	0	0	0	0	N1	N1	-1	0	0	N1	0	0
U7	0	0	-2N1	0	0	0	0	0	0	0	S1	0

**For Bottom Shim Plate**

INDEX (No.)	Display Shim Direction and Shim Location											
	0	30	60	90	120	150	180	210	240	270	300	330
L1	0	0	0	0	0	0	0	0	0	0	0	0
L2	0	0	0	0	0	0	0	0	0	0	0	0
L3	0	0	0	0	0	0	0	0	0	0	0	0
L4	0	0	0	0	0	0	0	0	0	0	0	0
L5	0	0	0	0	0	0	0	0	0	0	0	0
L6	0	0	S1	0	0	0	0	0	0	0	0	0
L7	N2	0	0	-3N1	0	0	N2	N1	N3	0	0	0


**HOW TO READ THE TABLE**  
ILLUSTRATION 29

Rev 4

### 6-4 Passive Shimming

#### Note

Refer to the following illustration to understand the position relationship between upper and lower shim plate.



**POSITION RELATIONSHIP OF UPPER AND LOWER SHIM PLATE  
ILLUSTRATION 30**

Rev 4

### 6-4-1 Upper Shim Plate Shim Setting

1. Remove upper shim plate from magnet by loosening A and B screws.


**Note**

You can remove the Shim Plates without removing the Mapping Fixture



**CAUTION**

When upper shim plates are installed to magnet, align 0degree of shim plates to left side(view from under side)as illustration. You can install center shim plate by rotating 90degree to magnet. The each outer shim plates can be installed to reverse position.




**UPPER SHIM PLATE 1**  
ILLUSTRATION 31

Rev 4

### 6-4-1 Upper Shim Plate Shim Setting (continued)

2. Remove cover from upper shim plate by loosening C screws.



**UPPER SHIM PLATE 2**  
ILLUSTRATION 32


Rev 4

**6-4-1 Upper Shim Plate Shim Setting (continued)**

3. Verify to check the pole(N-pole or S-pole) of the Shim by using N-S Checker.


**For checking N-axis:**

Attach N-S Checker to the Shim as shown. The magnet is repeling each other so the Shim facing the N-S Checker is N-Axis.



**NS CHECKER**  
ILLUSTRATION 33

4. Add or take out Shim Magnet on the Shim Board according to the output from the pshim program.



Cross Section of the Shim Board

**ADD SHIM MAGNET**  
ILLUSTRATION 34

5. Installtion is reverse order of removal.


Rev 4

### 6-4-2 Lower Shim Plate Shim Setting

1. Remove upper shim plate from magnet by loosening A and B screws.



When lower shim plates are install to magnet,align 0degree of shim plates to right side (view from upper side) as illustration. You can install center shim plate by rotating 90degree to magnet.The each outer shim plates can be installed to reverse position.




**LOWER SHIM PLATE 1**  
ILLUSTRATION 35

Rev 4

### 6-4-2 Lower Shim Plate Shim Setting (continued)

2. Remove cover from lower shim plate by loosening C screws.




Rev 4

### 6-4-2 Lower Shim Plate Shim Setting (continued)

3. Verify to check the pole(N-pole or S-pole) of the Shim by using N-S Checker.


**For checking N-axis:**

Attach N-S Checker to the Shim as shown. The magnet is repeling each other so the Shim facing the N-S Checker is N-Axis.



**NS CHECKER  
ILLUSTRATION 37**

4. Add or take out Shim Magnet on the Shim Board according to the output from the pshim program.



Cross Section of the Shim Board

**ADD SHIM MAGNET  
ILLUSTRATION 38**

5. Installtion is reverse order of removal.

Rev 4

**Revision History**

<b>Rev</b>	<b>Date</b>	<b>Author</b>	<b>Primary Reasons For Change</b>
0	Dec 20, 2000	Y. Masumo	Initial Version
1	May 19, 2001	Y. Masumo	Misc Correction
2	Oct 16, 2001	Y. Masumo	Misc Correction
3	Jan 28, 2002	Y. Masumo	Misc Correction
4	Jan 28, 2002	Y. Masumo	Page 4 to 9: Clarify the character and illustration. P39: Added to record Center Frequency Value. Page 42: Added explanation of upper and lower shim plate position relationship.