


# TABLE OF CONTENTS

TABLE OF CONTENTS .....	1
1- Phantom Setting .....	2
2- SCAN .....	2
3- Data Analysis .....	7
4- Data Sheet .....	10


## 1. Phantom Setting

1. Set the Head Coil on the Cradle.



**HEAD COIL POSITIONING**  
ILLUSTRATION 1


2. Set the DQA phantom and phantom holder according to the following illustration.



**DQA PHANTOM SETTING**  
ILLUSTRATION 1-1

## 2. SCAN


2-1 Flow Chart



**SCAN FLOWCHART**  
ILLUSTRATION 2-1

2-2 Patient Information

1. Select Patient Information and Protocol Icon.



Patient Information and Protocol Icon  
ILLUSTRATION 2-1

2. Click [New Pt] button in "PATIENT REGISTER".



[New Pt] button  
ILLUSTRATION 2-2

3. Input the following data on "patient information".

Patient Data	: geservice
Weight	: 50 (Kg)


The screenshot shows a form titled "PATIENT INFORMATION". The fields are as follows:
 

- Patient ID:
- Patient Name:
- Birth Date:  Age:  Sex:
- Weight:  (Lb)  (Kg)
- Rad:  Refer:
- Req Number:  Status:
- Description:
- History:
- Landmark:

Patient Information Input  
ILLUSTRATION 2-3

**2-3 Align/Send into the Magnet Center**

1. Turn alignment light ON.
2. Advance the cradle to the position where the axial A-light beam hits at 20mm outside from the axial Alignment line of Head Coil. Refer to the illustration.
3. Move the cradle to the position where the sagital A-light beam hits at the sagital Alignment line of Head Coil. Refer to the illustration



**ALIGNMENT LIGHT  
ILLUSTRATION 2**

4. Landmark in the sagittal and axial planes.
5. Press **[Adv to Scan]** button to send the phantom into the Magnet center.

**2-4 Patient Protocols**

1. Click **[Patient Position]**.
2. Enter the following Information.

**Table 1:**

Item	Value
Patient Position	Supine
Patient Entry	Head First
Coil	Head
Plane	Axial
Mode	2D
Pulse Seq	Spin Echo
Imaging Option	VBw + EDR

3. Enter the protocols according to the following list.

**Table 2:**

Item	Value
No of Echo	1
TE	25

**Table 2:**

Item	Value
TR	500
Band Width	10.42
FOV Size	25
Slice Thickness	5
Spacing	0
S/I Start	0
End	0
No of Slice	1
L/R Center	0
Freq	256
Phase	256
Nex	1
Phase FOV	1
Freq Dir	R/L
Auto Centrer Freq	Peak
AutoShim	ON

4. Select **Accept**.
5. Select **Save Series**.



Save Series  
ILLUSTRATION 2-1

**2-5 Auto Prescan**

1. Select **Auto Prescan** button



Auto Prescan

ILLUSTRATION 2-1


**2-6 Manual Prescan**

1. Select **Manual Prescan** button.



Manual Prescan  
ILLUSTRATION 2-1

2. Click on "SCAN TR" and check that the projection is displayed.



Projection  
ILLUSTRATION 2-2

**2-7 Scan**

1. Select **Scan** button



Scan  
ILLUSTRATION 2-1

2. Verify that scanning has been started.

**3. Data Analysis**


**3-1 Subtract Noise Image from Scan Image.**

1. Select Advantage Window Icon.



AW Icon  
ILLUSTRATION 3-1


2. Click [Add/Sub] button.



Add-Sub Button  
ILLUSTRATION 3-2

3. Subtract Image  
1st Image - 2nd Image

1. Click on "-".
2. Click on "Accept Negative
3. Select the Scan Image 1.
4. Click on "Select Set".(Left side)
5. Select the Scan Image 2.
6. Click "Select Set".(Right Side)
7. Click on "=".
8. Click "Quit".
9. The subtracted Image is saved as **Series 100\*\***



How to Subtract  
ILLUSTRATION 3-3


3-2 Measure the Mean and SD of 1st image.

1. Open the 1st Image by using the Miniviewer.
  - a. Select Advantage Window Icon.




AW Icon  
ILLUSTRATION 3-1

- b. Select "Exam", "Series", and "Image" from the patient list and then press [MiniViewer].




Miniviewer  
ILLUSTRATION 3-2

2. From "Cursors and Measuring", select ellipse and set the area to (3000±100mm<sup>2</sup>)



Elliptical Cursor  
ILLUSTRATION 3-3

3. Measure the Mean and SD.
  4. Record the Data into the Data Sheet.
  5. Press the keys at the same time to copy the current elliptical cursor :  
[Ctrl + Alt + F6] or [Ctrl + c].
- 3-3** Measure the Mean and SD of 2nd image.
1. Open the 2nd Image
    - a. Click on "Image +" button, and next image will come up.



Next Image  
ILLUSTRATION 3-1

2. Press the keys at the same time to paste the elliptical cursor:  
[Ctrl + Alt + F8] or [Ctrl + v].
  3. Measure the Mean and SD.
  4. Record the Data into the Data Sheet.
- 3-4** Measure the Subtracted Image.
1. Open the Subtracted Image.
    - a. Click on "Image +" button, and next image will come up.
  2. Measure the Mean and SD of the Subtracted Image.
  3. Record the Data into the Data Sheet.

**4. Data Sheet**

1. Record the Mean and SD values, and calculate the SNR:

$$SNR = (M1 \times \sqrt{2}) / SD(s)$$

where, M1 is the 1st Image Mean Value and SD(s) is the Subtracted Image SD value

TABLE 4-1

**MEAN AND SD VALUES**

	MEAN	SD
1st Image		
2nd Image		
Subtracted Image	-----	

TABLE 4-2

SNR VALUE

Calculated SNR	SPEC

I