

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- OVERVIEW	3
2- REQUIRED TOOLS	4
3- BODY AND HEAD BASELINE LOOPBACK DATA TESTS	6
3-1 Initial Setup Requirements	6
3-2 TPS Loopback Baseline Data Test	8
3-2-1 TPS Loopback Stability and Linearity Control Variables.....	11
3-3 Body Sense Loop Baseline Data Test.....	12
3-3-1 Body Sense Loop Test Stability and Linearity Control Variables	15
3-4 Head Sense Loop Baseline Data Test	16
3-4-1 Head Sense Loop Test Stability and Linearity Control Variables	20
3-5 Body Coil Load Baseline Data Test.....	21
3-5-1 Body Coil Load Test Stability and Linearity Control Variables.....	24
3-6 Body Dummy Load Baseline Data Test.....	25
3-6-1 Body Dummy Load Test Stability and Linearity Control Variables.....	28
3-7 Head Coil Load Baseline Data Test	29
3-7-1 Head Coil Load Test Stability and Linearity Control Variables	31
3-8 Head Dummy Load Baseline Data Test	32
3-8-1 Head Dummy Load Test Stability and Linearity Control Variables	34
4- COIL TUNING SENSE LOOP TESTS	35
4-1 Body Coil Tuning Sense Loop Test	35
4-2 Head Coil Tuning Sense Loop Test.....	41
5- BODY AND HEAD TROUBLESHOOTING TESTS	48
5-1 Body Dummy Load (Body Helix Output).....	48
5-2 Body Dummy Load (Body T/R Switch & Hybrid Output).....	51
5-3 Head Dummy Load (Head Helix Output).....	54
5-4 Head Dummy Load (Head T/R Output) (Head Helix Output)	57
5-5 Coil Tune Preamp Check.....	60
5-6 Preamp Bandwidth Check	63
6- RESTORATION CHECKLIST	66
7- RFT RESULTS	67
7-1 RFT Analysis	67
7-2 Report Program	67
7-3 RFT Header Data	70
7-4 RFT Sample Test Outputs.....	71
7-4-1 Low Power Magnitude (No Gradients) Stability Data.....	71
7-4-2 Low Power Phase (No Gradients) Stability Data	72
7-4-3 High Power Magnitude (No Gradients) Stability Data.....	73
7-4-4 High Power Phase (No Gradients) Stability Data	74
7-4-5 Low Power Magnitude (With Gradients) Stability Data.....	75
7-4-6 Low Power Phase (With Gradients) Stability Data	76
7-4-7 High Power Magnitude (With Gradients) Stability Data	77
7-4-8 High Power Phase (With Gradients) Stability Data	78
7-4-9 Linearity Data.....	79
7-4-10 Magnitude Droop Data.....	80
7-4-11 Sinc 2 Pulse Fidelity Data	81
7-4-12 Sinc 8 Pulse Fidelity Data	82
7-4-13 Bandwidth Data	83
7-4-14 Coil Tune Check Data	84
7-4-15 Config File Data	85
8- RFT THEORY	86
8-1 RFT Overview.....	86
8-2 Data Collection	87
8-3 Magnitude Stability Analysis (groups 1, 3, 5 and 7).....	89
8-4 Phase Stability Analysis (groups 2, 4, 6 & 8).....	89
8-5 Linearity And Carrier Leakage Analysis (group 9).....	90

8-6 Magnitude Droop Analysis (group 10) 91
8-7 Sinc 2 Pulse Fidelity Analysis (group 11) (TPS Loopback only) 91
8-8 Sinc 8 Pulse Fidelity Analysis (group 12) (TPS Loopback only) 91
8-9 Magnitude Bandwidth Analysis (group 13 & 14) 92
APPENDIX A - DATA SHEET 94
APPENDIX B - CV DESCRIPTIONS 95
REVISION HISTORY 97

1- OVERVIEW

RFT is a RF loopback test that quickly evaluates stability, magnitude and phase linearity, droop, pulse fidelity, and bandwidth of the RF transmit chain in a Signa system.

It is recommended that TPS Loopback Baseline Data Test is performed first. Run Body Sense Loop Baseline Data Test and Head Sense Loop Baseline Data Test next.

The receive line is not verified in this procedure.

For additional information, refer to Section 7-1 RFT Analysis or 8-1 RFT Overview.

IMPORTANT!

The RFT linearity specification is different for each setting of the EFB Bypass CV. A change will be implemented (via a patch from the ONLINE Center for 1.0T and 1.5T pre 8.3 software) so that the RF Test will use the RF Amplifier Type configuration file selection (Tube type or a Solid State type). The analysis portion of RFT will then enable the tighter linearity specification which is required when the configuration file is Solid State RF Amp Type.

Note

Disable the TNF (Transient Noise Filter) unit before running any RFT test. The TNF is located inside the rear door of the System Cabinet. Remember to enable the unit after RFT testing has completed

2- REQUIRED TOOLS

Refer to Tables 2-1, 2-2, 2-3, and 2-4 for the kits required and their contents.

TABLE 2-1
UNIVERSAL SST KIT 46-320383G1 — G7 CONTENTS

Item	Description	Part Number	Quantity
1	Cable, 5ft., RG58, with BNC Connectors	46-251710G4	5
2	Cable, 1ft., RG58, with BNC Connectors	46-282803G16	2
3	Cable, 0.5 ft., RG58, with BNC Connectors	46-251920G4	1
4	Cable, 5ft., RG188, with BNC Male/SMB Female	46-301549P5	1
5	DC Power Supply for Switchable Attenuator	46-320385P1	1
6	Switchable Attenuator	46-288458G1	1
7	1.0T RF Sense Loop (Only in G4 & G6 Kits)	46-321434G1	1
8	1.5T RF Sense Loop (Only in G1/G2/G3/G4/G5 Kits)	46-288446G1	1
9	Head Carriage Sense Loop "L" Bracket	2102499	1
10	Attenuator, Variable 0-10 dB (1dB steps)	46-255838P3	1
11	Attenuator, Variable 0-1 dB (0.1dB steps)	46-255838P7	1
12	Velcro Pads	46-307152P1	3
13	Case of Connectors	46-301042G1	1

TABLE 2-2
RF POWER MEASUREMENT KIT CONTENTS

Item	Description	Part Number
1.	RF Power Measurement Kit NOTE: G1 kit does not contain the 30 dB Load. 50 ohm, 200 Watt, 30dB Attenuator Bird Model 8322	46-317724G1 or G2 46-255837P10 or 46-317724P14
2.	RF Cable Test Kit NOTE: For 1.5T RF/PDU Cabinet: N Female to N Female, 50 ohm Type, Adapter	46-255816G1 46-265875P2
3.	0-10 dB Rotary Step Attenuator, (variable 1 dB steps)	46-255838P1

TABLE 2-3
**QTUNE KIT (46-287308G1/G2/G3) CONTENTS
 (USED FOR RFT COIL TUNING CHECK ONLY)**

Item	Description	Part Number	Quantity
1	Sense Coil Holder (for Head Coil) (G1 Kit only)	46-287006P1	1
2	Sense Coil Holder (for Head Coil) (G2 & G3 Kit only)	46-307899P1	1
3	1.0T RF Sense Loop (Only in G3 Kit)	46-321434G1	2
4	1.5T RF Sense Loop (Only in G1/G2/G3 Kits)	46-288446G1	2

TABLE 2-4
OTHER ITEMS NEEDED

Item	Description	Part Number	Quantity
1	50-ohm, 30 dB, 200 Watt, Dummy Load NOTE: Dummy Load (Body T/R Switch & Hybrid Output) test requires 2 (two) Dummy Loads.	46-255837P10	1
2	Grafidy Holder Base Plate	46-271410G1	1
3	Attenuator, Variable 0-70 dB (10dB steps)	46-255838P2	1
4	N-Shorting Caps *This part is supplied with the Coil Tune Extension kit	* 46-265916P1	2

3- BODY AND HEAD BASELINE LOOPBACK DATA TESTS

3-1 Initial Setup Requirements




EQUIPMENT DAMAGE POSSIBILITY! The Switchable Attenuator contains ferrous switch relays. Do not take the attenuator into the magnet room; doing so can cause damage to the switch relays.

1. Install the Switchable Attenuator, 46-288458G1, on the rear of the Systems Cabinet as shown in Illustration 3-1.




EQUIPMENT DAMAGE! CONNECT SWITCHABLE ATTENUATOR BEFORE CONNECTING SENSE LOOP TO PREVENT DAMAGE TO TPS RECEIVER BOARD. ALSO, LEAVE SWITCH S1 (PREAMP BIAS) ON SWITCHABLE ATTENUATOR IN "OUT" (DEFAULT) POSITION OR THE INTERNAL PROTECTIVE DC BLOCKER IS BYPASSED. (WITH S1 OUT, DS1 "BIAS ON" LED SHOULD BE OFF AFTER SWITCHABLE ATTENUATOR IS CONNECTED.



SWITCHABLE ATTENUATOR
ILLUSTRATION 3-1

2. On some Switchable Attenuators: Push in S1, then remove J5 power to restore the DS1 to OFF (DS1 latches).

3. Disable the TNF (Transient Noise Filter) unit before running any RFT test. The TNF is located inside the rear door of the System Cabinet.
4. See Illustration 3-2. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.




SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
ILLUSTRATION 3-2


3-2 TPS Loopback Baseline Data Test

This test verifies the Exciter and Receiver hardware is good. TPS Loopback Baseline Data Test is performed in the Body mode. Results from this test are referenced automatically by subsequent RFT tests, therefore, this test is performed first. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-2-1, Table 3-2.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. Install the Switchable Attenuator, 46-288458G1, at the rear of the Systems Cabinet as shown in Illustration 3-3A. **Or** install the 0 to 10 dB Variable Attenuator, 46-255838P3 and the DC Block, 46-301549P15, at the rear of the Systems Cabinet referring to Illustration 3-3B (Alternate Hardware Setup). This alternate test setup will perform the TPS Loopback test without the use of the Switchable Attenuator. The alternate hardware setup is useful if attempting to identify a hardware problem associated with the Switchable Attenuator.



TPS LOOPBACK SETUP
 ILLUSTRATION 3-3A



ALTERNATE HADWARE SETUP FOR TPS LOOPBACK
 ILLUSTRATION 3-3B

3. Prepare the system to scan in Body mode per Table 3-1.

TABLE 3-1
BODY SCAN PROTOCOL: TPS LOOPBACK BASELINE DATA TEST

—Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen].

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **0**<ENTER>.

CV Name: **Sense Loop**, CV Value: **0**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **1**<ENTER>.

CV Name: **BW Cal**, CV Value: **1**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 15

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

4. **[Manual Prescan].**

5. From the “Windows” Menu bar on the Manual Prescan window, select “**Two Windows**”.


Plot Type = **Magnitude**

Plot Gain = **1**

Plot Type = **Magnitude**

Plot Gain = **1**

There are 4 Frames of rectangular pulses. Two of these pulses will be viewed, one in each window.



Note: This basic signal will be displayed during any of the tests on the following pages while in prescan.

Note: If the signal is low (~2%), check the AC Adapter Power Supply for power to the Switchable Attenuator.

L17269

**TPS LOOPBACK MANUAL PRESCAN SIGNAL DISPLAY
ILLUSTRATION 3-4**

6. Set TG to 200.

IMPORTANT!

For TPS Loopback Set Up ONLY: See Illustration 3-4. The R2 Power Spectrum signal should be between 40% & 90%. This may require attenuator A to be set at 0 dB.

7. Select **[Done]**.

8. Select **[Scan]**. Scan time is approximately 5 minutes and 30 seconds.

9. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in the appendix.

10. If not proceeding to next sub-section, perform System Restoration.

3-2-1 TPS Loopback Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-2:


**TABLE 3-2
 TPS LOOPBACK CONTROL VARIABLES**

User CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFampMode	0	0	0
Sense Loop	0	0	0
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
Stab/Grad	4	4	0
Bandwidth (see *)	1	0	1
BW Cal (see *)	1	0	0


3-3 Body Sense Loop Baseline Data Test

This test evaluates the health of the Body Transmit chain and the RF field in the bore. If this test passes the Body Transmit chain is good. TPS Loopback Baseline Data Test is performed first. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-3-1, Table 3-4.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. Refer to Illustrations 3-5 and 3-6 for the hardware setup.



BODY SENSE LOOP SET-UP
 ILLUSTRATION 3-5



BODY COIL HARDWARE SET-UP
 ILLUSTRATION 3-6

3. Prepare the system to scan in Body mode per Table 3-3.

TABLE 3-3
BODY SCAN PROTOCOL: BODY SENSE LOOP BASELINE

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark per Illustration, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen].

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **4**<ENTER>.

CV Name: **Sense Loop**, CV Value: **1**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **1**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

4. **[Manual Prescan]**.

5. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".

Plot Type = **Magnitude**

Plot Gain = **1**

Plot Type = **Magnitude**

Plot Gain = **1**

6. Advance TG to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
7. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
8. **[Done]**.
9. **[Scan]**.
10. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in the appendix.
11. If not proceeding to next subsystem, perform System Restoration.

3-3-1 Body Sense Loop Test Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-4:

TABLE 3-4
BODY SENSE LOOP CONTROL VARIABLES

User CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFAmpMode	4	4	4
Sense Loop	1	1	1
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
Stab/Grad	4	4	0
Bandwidth (see*)	1	0	1
BW Cal	0	0	0

3-4 Head Sense Loop Baseline Data Test

This test evaluates the health of the Head Transmit chain and the RF field in the bore. If this test passes the Head Transmit chain is good. TPS Loopback Baseline Data Test is performed first. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-4-1, Table 3-6.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. Prepare the system to scan in Head mode per Table 3-5.

TABLE 3-5
HEAD SCAN PROTOCOL: HEAD SENSE LOOP BASELINE DATA TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark per Illustration, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.2**<ENTER>(o=Other, 16.2 =series) to load the head protocol

OR select **other** and select protocol **16** (RF test)and select series **2** (Head scan).

C. In the Additional Parameters field:

[User CVs Screen] .

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **2**<ENTER>.

CV Name: **Sense Loop**, CV Value: **1**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **0**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>


+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

- 3. See Illustration 3-7, 3-8 and 3-9 for the head setup.
- A. Attach Velcro as shown in Illustration 3-7.




VELCRO ATTACHEMENT
ILLUSTRATION 3-7

- B. Position sense loop per Illustration 3-8. DO NOT use tuning ring. Position the sense loop on the sense loop bracket as shown.


WARNING!

FOR THIS TEST CONFIGURATION, DO NOT PLACE THE SENSE LOOP IN THE HEAD COIL. THE HEAD COIL ASSEMBLY MAY BE DAMAGED.



POSITIONING SENSE LOOP AT HEAD COIL
ILLUSTRATION 3-8

C. Position and landmark Head TLT Sphere in the Loader per Illustration 3-9. DO NOT use tuning ring.



POSITIONING HEAD TLT SPHERE INSIDE THE HEAD LOADER
ILLUSTRATION 3-9

3-4-1 Head Sense Loop Test Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-6:


TABLE 3-6
HEAD SENSE LOOP CONTROL VARIABLES

User CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFampMode	2	2	2
Sense Loop	1	1	1
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
Stab/Grad	4	4	0
Bandwidth	0	0	0
BW Cal	0	0	0

3-5 Body Coil Load Baseline Data Test

This test evaluates the health of the Body Transmit chain out to the Body Coil with the RF Amplifier loaded by the coil, however, it does not evaluate the RF field in the bore. TPS Loopback Baseline Data Test is performed first. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-5-1, Table 3-8.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. See Illustration 3-11 for the hardware setup.



BODY COIL LOAD LOOPBACK SETUP
 ILLUSTRATION 3-11

3. Prepare the system to scan in Body mode per Table 3-7.

TABLE 3-7
BODY SCAN PROTOCOL: BODY COIL LOAD BASELINE DATA TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen].

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **4**<ENTER>.

CV Name: **Sense Loop**, CV Value: **0**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **0**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

4. **[Manual Prescan]**.

5. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".

Plot Type = **Magnitude**

Plot Gain = **1**

Plot Type = **Magnitude**

Plot Gain = **1**

6. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 over-range. Do NOT adjust R1 or R2.
7. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
8. Select [**Done**].
9. Select [**Scan**].
10. When scan completes, view results by using the Report Manager tool. Record data for each test Data Sheet 0.7T found in Appendix A.
11. If not proceeding to next sub-section, perform System Restoration.

3-5-1 Body Coil Load Test Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-8:


**TABLE 3-8
 BODY COIL LOAD CONTROL VARIABLES**

CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFampMode	4	4	4
SenseLoop	0	0	0
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
StabGrad	4	4	0
Bandwidth	0	0	0
BW Cal	0	0	0

3-6 Body Dummy Load Baseline Data Test

This test evaluates the Body Transmit chain up to the Body RF Out at the RF/PDU cabinet with the RF Amplifier loaded by a 50 ohm Dummy Load. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-6-1, Table 3-10.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. See Illustration 3-12 for the hardware setup.



BODY DUMMY LOAD LOOPBACK SETUP
 ILLUSTRATION 3-12

3. Prepare the system to scan in Body mode per Table 3-9.

TABLE 3-9
BODY SCAN PROTOCOL: BODY DUMMY LOAD BASELINE

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen]

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **3**<ENTER>.

CV Name: **Sense Loop**, CV Value: **0**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **0**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

4. **[Manual Prescan].**
5. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".
Plot Type = **Magnitude**
Plot Gain = **1**
Plot Type = **Magnitude**
Plot Gain = **1**
6. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
7. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
8. **[Done].**
9. **[Scan].**
10. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in the appendix.
11. If not proceeding to next sub-section, perform System Restoration.

3-6-1 Body Dummy Load Test Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-10:

TABLE 3-10
BODY DUMMY LOAD CONTROL VARIABLES

User CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFAmpMode	3	3	3
Sense Loop	0	0	0
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
Stab/Grad	4	4	0
Bandwidth	0	0	0
BW Cal	0	0	0

3-7 Head Coil Load Baseline Data Test

This test evaluates the health of the Head Transmit chain out to the Head Coil with the RF Amplifier loaded by the coil, however, it does not evaluate the RF field in the bore. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-7-1, Table 3-12.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. Prepare the system to scan in Head mode per Table 3-11.

TABLE 3-11
HEAD SCAN PROTOCOL: HEAD COIL LOAD BASELINE DATA TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, , landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.2**<ENTER>(o=Other, 16.2 =series) to load the head protocol

OR select **other** and select protocol **16** (RF test) and select series **2** (Head scan).

C. In the Additional Parameters field:

[User CVs Screen]

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **2**<ENTER>.

CV Name: **Sense Loop**, CV Value: **0**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidly**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **0**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations] [Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>


+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations] [Download]**.

3. See Illustration 3-13 for the hardware setup.



HEAD COIL LOAD LOOPBACK SETUP
ILLUSTRATION 3-13

4. **[Manual Prescan].**

5. From the “Windows” Menu bar on the Manual Prescan window, select “**Two Windows**”.

Plot Type = **Magnitude**

Plot Gain = **1**

Plot Type = **Magnitude**

Plot Gain = **1**

6. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 over-range. Do NOT adjust R1 or R2.

7. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.

8. Select **[Done]**.

9. Select **[Scan]**.

10. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in Appendix A.

11. If not proceeding to next sub-section, perform System Restoration.

3-7-1 Head Coil Load Test Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-12:

TABLE 3-12
HEAD COIL LOAD CONTROL VARIABLES

User CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFAmpMode	2	2	2
Sense Loop	0	0	0
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
Stab/Grad	4	4	0
Bandwidth	0	0	0
BW Cal	0	0	0

3-8 Head Dummy Load Baseline Data Test

This test evaluates the Head Transmit chain up to the Head RF Out at the RF/PDU cabinet with the RF Amplifier loaded by a 50 ohm Dummy Load. Control Variables for Stability and Linearity trouble-shooting are listed in Sub-Section 3-8-1, Table 3-14.

1. Prepare the system per sub-section 3-1 (TNF and SSM).
2. Prepare the system to scan in Head mode per Table 3-13.

TABLE 3-13

HEAD SCAN PROTOCOL: HEAD DUMMY LOAD BASELINE DATA TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, , landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.2**<ENTER>(o=Other, 16.2 =series) to load the head protocol

OR select **other** and select protocol **16** (RF test) and select series **2** (Head scan).

C. In the Additional Parameters field:

[User CVs Screen]

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **1**<ENTER>.

CV Name: **Sense Loop**, CV Value: **0**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **1**<ENTER>.

CV Name: **Lin&Fidly**, CV Value: **1**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **4**<ENTER>.

CV Name: **Bandwidth**, CV Value: **0**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations] [Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>


+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations] [Download]**.

3. See Illustration 3-14 for the hardware setup.



HEAD MODE, DUMMY LOAD LOOPBACK SETUP
ILLUSTRATION 3-14

4. **[Manual Prescan].**

5. From the “Windows” Menu bar on the Manual Prescan window, select “**Two Windows**”.

Plot Type = **Magnitude**

Plot Gain = **1**

Plot Type = **Magnitude**

Plot Gain = **1**

6. Advance TG to 200, while adjusting variable attenuators A and B to prevent R1 and R2 over-range. Do NOT adjust R1 or R2.

7. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.

8. Select **[Done]**.

9. Select **[Scan]**.

10. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in Appendix A.

11. If not proceeding to next sub-section, perform System Restoration.

3-8-1 Head Dummy Load Test Stability and Linearity Control Variables

1. Control Variables for Stability and Linearity trouble-shooting are listed in Table 3-14:

TABLE 3-14
HEAD DUMMY LOAD CONTROL VARIABLES


User CV Name	Baseline	Stability T/S	Linearity T/S
EFB Bypass	0	0	0
RFAmpMode	1	1	1
Sense Loop	0	0	0
StabNoGrad	1	1	0
Lin&Fidly	1	0	1
Stab/Grad	4	4	0
Bandwidth	0	0	0
BW Cal	0	0	0

4- COIL TUNING SENSE LOOP TESTS

The Coil Tuning Sense Loop Tests evaluate the coils (body or head) resonant frequency and quality factor ($q=cf/bw$).

4-1 Body Coil Tuning Sense Loop Test

1. Setup the initial hardware using the following steps:
 - A. Bypass the TNF unit at the rear of the System Cabinet.
 - B. Refer to Illustration 4-1, at the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.




SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
ILLUSTRATION 4-1

2. For an accurate Body Coil Tuning reading, the Body Hybrid Splitter must be electrically isolated from the Body Coil. Remove the blue I and Q RF lines from J3 and J4 of the Body Hybrid Splitter.

- A. Once the I and Q lines have been removed, attach an N Shorting Cap (46-265916P1) to each line. Refer to Illustration 4-2.


Note

The N Shorting cap is supplied with the Coil Tune Extension kit (46-265917G1).



TERMINATING I & Q RF LINES WITH SHORTING CAPS
 ILLUSTRATION 4-2

- 3. Connect system as shown in Illustration 4-3 and 4-4.



BODY COIL TUNE SENSE LOOP CHECK
 ILLUSTRATION 4-3


- A. Position the Sense coils per Illustration 4-4 using sticky backed velcro.
- B. RF Coil type identified by part number (not description) on the sticker at rear of the Coil Assembly, see Table 4 for help:

TABLE 4
SENSE COIL POSITIONING

RF Coil Type	Part Number	Measure in from RF Coil Tube edge:	Measured Value
BRM or BRM/D	46-328500G1 46-328500G2	Front edge of RF Coil Tube (white). Cone edge and RF Coil Tube meeting place.	In approximately 16 inches (approximately 408 mm)
CRM	2186651	Front edge of RF Coil Tube (white). Cone edge and RF Coil Tube meeting place.	In approximately 10.5 inches (approximately 270 mm)
BRM/F	2218111 2218111-2	Front edge of RF Coil Tube (white). Cone edge and RF Coil Tube meeting place.	In approximately 0.9 inches (in tenths) or 7/8 inches (approximately 2.2 mm)

Note

When routing the cables for the Sense Coils, try to lay them as straight as possible and without crossing them. Crossed cables will affect data collection due to crosstalk. The sense coils must be positioned correctly or tuning problems will occur. Placing the sense coils too far into the Body Coil will cause over-coupling (creating 2 false peaks). Placing the sense coils too far out of the Body Coil will cause a decrease in frequency and sensitivity. Ensure that the sense coil foil run is placed the proper distance in per Table 4-4, Sense Coil Positioning).



SENSE COIL POSITIONING
 ILLUSTRATION 4-4

4. Prepare the system to scan in Body mode per Table 4-1.

TABLE 4-1
BODY SCAN PROTOCOL: BODY COIL TUNING SENSE LOOP TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen].

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **7**<ENTER>.

CV Name: **Sense Loop**, CV Value: **1**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **0**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **0**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **0**<ENTER>.

CV Name: **Bandwidth**, CV Value: **1**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>


+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Start **[Manual Prescan]**.
6. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".
Plot Type = **Magnitude**
Plot Gain =1
Plot Type = **Magnitude**
Plot Gain =1
7. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
8. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
9. Select **[Done]**.
10. Select **[Scan]**. When scan completes, view results by using the Report Manager tool. Record the "Center Frequency" value for the sense coil plane. (This value will later be used in calculating Peak Frequency).
11. Perform a new Body Coil Tune Sense Loop on all remaining planes. Refer to Illustration 4-5 for plane references (positions 2 through 4).



Ⓝ SENSE COIL (THE NUMBER INSIDE THE CIRCLE REPRESENTS THE TUNING POSITION, 1 - 4)

L18420

BODY SENSE COIL PLACEMENT
ILLUSTRATION 4-5

12. After obtaining Head Coil Tune Check Scans from all planes (0, 45, 90, and 315 degrees), view results by using the Report Manager tool. Record data for each test Data Sheet 0.7T found in the appendix.

- A. Select **Report Tool** under **[Utilities]** on the Service Desktop. Select the **File** menu, then select **Open** from the pull-down sub-menu. Choose the pertinent file from the *Load Data File* pop up window. The most recent files appear at the end of the file directory.
- B. **Center Frequency:** Record the center frequency obtained at each plane. Subtract the lowest frequency value from the highest frequency. The “mode delta” difference must be less than or equal to **150,000 Hz (1.5T) or 100,000 Hz (1.0T)**.

Example:

- degree plane: body coil tune check scan: 42.731 MHz
- degree plane: body coil tune check scan: 42.683 MHz
- $42.731\text{MHz} - 42.683\text{ MHz} = 48,000\text{ Hz}$ or 48 KHz. These frequency values are within spec.

13. If the bandwidth is not correct, then perform “Preamplifier Bandwidth Check” and troubleshoot the problem.

4-2 Head Coil Tuning Sense Loop Test

1. Setup the initial hardware using the following steps:


- A. Disconnect the 2 (two) white BNC cables and the Head Coil Quick Disconnect Adapter from the Head Coil.
- B. Locate and move both JP1 and JP2 jumpers to the Tune position (specific to the coil design). Refer to Illustrations 4-6 and 4-7.

Note

The 1.5T Isolation Boards and Hybrid Assembly are not removed, however, the semi-rigid coax's are removed.

WARNING!


FOR THIS TEST CONFIGURATION, DO NOT LEAVE THE HEAD COIL CONNECTED TO THE CARRIAGE COVER ASSEMBLY. THE HEAD T/R ASSEMBLY MAY BE DAMAGED. REMOVE THE BNC CABLES AND HEAD COIL QUICK DISCONNECT ADAPTER BEFORE INSTALLING THE HEAD COIL IN THE BORE.



JUMPER POSITIONS ON 1.5T HEAD COIL
ILLUSTRATION 4-6

JUMPER POSITIONS ON 1.0T HEAD COIL
ILLUSTRATION 4-7

2. See Illustration 4-8. On the front of the SSM place the DRV FLT switches as follows:
- TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
ILLUSTRATION 4-8

3. Bypass the TNF unit at the rear of the System Cabinet.

4. Prepare the system to scan in Head mode per Table 4-2.

TABLE 4-2
HEAD SCAN PROTOCOL: HEAD COIL TUNE SENSE LOOP CHECK TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, , landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.2**<ENTER>(o=Other, 16.2 =series) to load the head protocol

OR select **other** and select protocol **16** (RF test) and select series **2** (Head scan).

C. In the Additional Parameters field:

[User CVs Screen]

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **7**<ENTER>.

CV Name: **Sense Loop**, CV Value: **1**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **0**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **0**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **0**<ENTER>.

CV Name: **Bandwidth**, CV Value: **1**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>


+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Setup system per Illustration 4-9 and Illustration 4-10.




HEAD COIL TUNE SENSE LOOP CHECK
 ILLUSTRATION 4-9

6. Position the Sense coils per Illustration 4-10.
 - A. Place Head foam piece into front of head coil.
 - B. Place Sense Coils into foam positioner until the edge of the sense coil is aligned with the front edge of the Head Coil.
 - C. Landmark per Illustration 4-10.

Note


When routing the cables for the Sense Coils, try to lay them as straight as possible and without crossing them. Crossed cables will affect data collection due to crosstalk.



L1781A

HEAD COIL SENSE COIL POSITIONING
ILLUSTRATION 4-10

7. Start **[Manual Prescan]**.
8. From the "Windows" Menu bar on the Manual Prescan window, select **"Two Windows"**.
Plot Type = **Magnitude**
Plot Gain = **1**
Plot Type = **Magnitude**
Plot Gain = **1**
9. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
10. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
11. Select **[Done]**.
12. Select **[Scan]**.
13. For scan #1, place sense coils at position 1. For subsequent scans, place sense coils at next location. Refer to Illustration 4-11.



SENSE COIL POSITIONING / PLANE REFERENCES
ILLUSTRATION 4-11

14. Perform a new Head Coil Tune Sense Loop Test scan on all remaining planes (positions 2 through 4). Refer to Illustration 4-11 for plane / position references.

15. After obtaining Head Coil Tune Check Scans from all planes (0, 45, 90, and 315 degrees), view results by using the Report Manager tool. Record data for each test Data Sheet 0.7T found in the appendix.

- A. Select **Report Tool** under **[Utilities]** on the Service Desktop. Select the **File** menu, then select **Open** from the pull-down sub-menu. Choose the pertinent file from the *Load Data File* pop up window. The most recent files appear at the end of the file directory.
- B. **Center Frequency:** Record the center frequency obtained at each plane. Subtract the lowest frequency value from the highest frequency. The “mode delta” difference must be less than or equal to **150,000 Hz (1.5T) or 100,000 Hz (1.0T).**

Example:

- degree plane: body coil tune check scan: 42.731 MHz
- degree plane: body coil tune check scan: 42.683 MHz
- $2.731\text{MHz} - 42.683\text{ MHz} = 48,000\text{ Hz}$ or 48 kHz. These frequency values are within spec.


16. If the bandwidth is not correct, then perform “Pre-amplifier Bandwidth Check” and troubleshoot the problem.

5- BODY AND HEAD TROUBLESHOOTING TESTS

5-1 Body Dummy Load (Body Heliax Output)


1. Refer to Illustration 5-1. On the front of the SSM place the DRV FLT switches as follows:

- TR switch to the DIS (disable faults) position.
- DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
 ILLUSTRATION 5-1

2. Bypass the TNF unit at the rear of the System Cabinet.
3. Refer to Illustration 5-2 for hardware setup:



BODY MODE DUMMY LOAD LOOPBACK (BODY HELIAX OUTPUT)
 ILLUSTRATION 5-2

4. Prepare the system to scan in Body mode per Table 5-1.

TABLE 5-1
BODY SCAN PROTOCOL: BODY DUMMY LOAD (BODY HELIAX OUTPUT) TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen] (see Table 5-2 for CV Values):

CV Name: **EFB Bypass**, CV Value: <ENTER>.

CV Name: **RFampMode**, CV Value: <ENTER>.

CV Name: **Sense Loop**, CV Value: <ENTER>.

CV Name: **StabNoGrad**, CV Value: <ENTER>.

CV Name: **Lin&Fidlty**, CV Value: <ENTER>.

CV Name: **Stab/Grad**, CV Value: <ENTER>.

CV Name: **Bandwidth**, CV Value: <ENTER>.

CV Name: **BW Cal**, CV Value: <ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Control Variables for Stability and Linearity trouble-shooting are listed in Table 5-2 for **[User CV's]** screen select (refer to Appendix B for CV descriptions):


TABLE 5-2
BODY DUMMY LOAD (BODY HELIAX OUTPUT) CONTROL VARIABLES

User CV Name	Stability T/S	Linearity T/S
EFB Bypass	0	0
RFAmpMode	3	3
Sense Loop	0	0
StabNoGrad	1	0
Lin&Fidlty	0	1
Stab/Grad	4	0
Bandwidth	0	0
BW Cal	0	0

6. Start **[Manual Prescan]**.
7. From the "Windows" Menu bar on the Manual Prescan window, select **"Two Windows"**.
 Plot Type = **Magnitude**
 Plot Gain = **1**
 Plot Type = **Magnitude**
 Plot Gain = **1**
8. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
9. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
10. Select **[Done]**.
11. Select **[Scan]**. When scan completes, view results by using the Report Manager tool (Release 8.2) or the system PC FTP tool (Release 8.1). Record data for each test in Data Sheet 0.7T found in the appendix.


5-2 Body Dummy Load (Body T/R Switch & Hybrid Output)

1. Refer to Illustration 5-3. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
 ILLUSTRATION 5-3

2. Bypass the TNF unit at the rear of the System Cabinet.
3. Refer to Illustration 5-4 for hardware setup.



BODY MODE DUMMY LOAD LOOPBACK (BODY T/R & HYBRID OUTPUT) SETUP
 ILLUSTRATION 5-4

4. Prepare the system to scan in Body mode per Table 5-3.

TABLE 5-3

BODY SCAN PROTOCOL: BODY DUMMY LOAD (BODY T/R SWITCH & HYBRID OUTPUT) TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen] (see Table 5-2 for CV Values):

CV Name: **EFB Bypass**, CV Value: <ENTER>.

CV Name: **RFampMode**, CV Value: <ENTER>.

CV Name: **Sense Loop**, CV Value: <ENTER>.

CV Name: **StabNoGrad**, CV Value: <ENTER>.

CV Name: **Lin&Fidlty**, CV Value: <ENTER>.

CV Name: **Stab/Grad**, CV Value: <ENTER>.

CV Name: **Bandwidth**, CV Value: <ENTER>.

CV Name: **BW Cal**, CV Value: <ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Control Variables for Stability and Linearity trouble-shooting are listed in Table 5-4 for **[User CV's]** screen select (refer to Appendix B for CV descriptions):


TABLE 5-4
BODY DUMMY LOAD (BODY T/R SWITCH & HYBRID OUTPUT) CONTROL VARIABLES

User CV Name	Stability T/S	Linearity T/S
EFB Bypass	0	0
RFAmpMode	3	3
Sense Loop	0	0
StabNoGrad	1	0
Lin&Fidlty	0	1
Stab/Grad	4	0
Bandwidth	0	0
BW Cal	0	0

6. Start **[Manual Prescan]**.
7. From the "Windows" Menu bar on the Manual Prescan window, select **"Two Windows"**.
 Plot Type = **Magnitude**
 Plot Gain = **1**
 Plot Type = **Magnitude**
 Plot Gain = **1**
8. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
9. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
10. Select **[Done]**.
11. Select **[Scan]**. When scan completes, view results by using the Report Manager tool. Record data for each test Data Sheet 0.7T found in the appendix.


5-3 Head Dummy Load (Head Heliax Output)

1. Refer to Illustration 5-5. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
 ILLUSTRATION 5-5

2. Bypass the TNF unit at the rear of the System Cabinet.
3. Refer to Illustration 5-6 for hardware setup.



HEAD MODE DUMMY LOAD LOOPBACK (HEAD HELIAX OUTPUT)
 ILLUSTRATION 5-6

4. Prepare the system to scan in Head mode per Table 5-5.

TABLE 5-5
HEAD SCAN PROTOCOL: HEAD DUMMY LOAD (HEAD HELIAX OUTPUT) TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.2**<ENTER>(o=Other, 16.2 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **2** (RFT Head).

C. In the Additional Parameters field:

[User CVs Screen] (see Table 5-2 for CV Values):

CV Name: **EFB Bypass**, CV Value: <ENTER>.

CV Name: **RFampMode**, CV Value: <ENTER>.

CV Name: **Sense Loop**, CV Value: <ENTER>.

CV Name: **StabNoGrad**, CV Value: <ENTER>.

CV Name: **Lin&Fidlty**, CV Value: <ENTER>.

CV Name: **Stab/Grad**, CV Value: <ENTER>.

CV Name: **Bandwidth**, CV Value: <ENTER>.

CV Name: **BW Cal**, CV Value: <ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Control Variables for Stability and Linearity trouble-shooting are listed in Table 5-6 for **[User CV's]** screen select (refer to Appendix B for CV descriptions):


TABLE 5-6
 HEAD DUMMY LOAD (HEAD HELIAX OUTPUT) CONTROL VARIABLES

User CV Name	Stability T/S	Linearity T/S
EFB Bypass	0	0
RFAmpMode	1	1
Sense Loop	0	0
StabNoGrad	1	0
Lin&Fidlty	0	1
Stab/Grad	4	0
Bandwidth	0	0
BW Cal	0	0

6. Start **[Manual Prescan]**.
7. From the "Windows" Menu bar on the Manual Prescan window, select **"Two Windows"**.
 Plot Type = **Magnitude**
 Plot Gain = **1**
 Plot Type = **Magnitude**
 Plot Gain = **1**
8. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
9. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
10. Select **[Done]**.
11. Select **[Scan]**. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in Appendix A.


5-4 Head Dummy Load (Head T/R Output) (Head Heliax Output)

1. Refer to Illustration 5-7. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
 ILLUSTRATION 5-7

2. Bypass the TNF unit at the rear of the System Cabinet.
3. Refer to Illustration 5-8 for the hardware setup:



HEAD MODE DUMMY LOAD LOOPBACK (HEAD T/R) SETUP
 illustration 5-8

4. Prepare the system to scan in Head mode per Table 5-7.

TABLE 5-7

HEAD SCAN PROTOCOL: HEAD DUMMY LOAD (HEAD T/R OUTPUT) (HEAD HELIAX OUTPUT) TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.2**<ENTER>(o=Other, 16.2 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **2** (RFT Head).

C. In the Additional Parameters field:

[User CVs Screen] (see Table 5-2 for CV Values):

CV Name: **EFB Bypass**, CV Value: <ENTER>.

CV Name: **RFampMode**, CV Value: <ENTER>.

CV Name: **Sense Loop**, CV Value: <ENTER>.

CV Name: **StabNoGrad**, CV Value: <ENTER>.

CV Name: **Lin&Fidlty**, CV Value: <ENTER>.

CV Name: **Stab/Grad**, CV Value: <ENTER>.

CV Name: **Bandwidth**, CV Value: <ENTER>.

CV Name: **BW Cal**, CV Value: <ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Control Variables for Stability and Linearity trouble-shooting are listed in Table 5-8 for **[User CV's]** screen select (refer to Appendix B for CV descriptions):


TABLE 5-8
 HEAD DUMMY LOAD (HEAD T/R OUTPUT) (HEAD HELIAX OUTPUT) CONTROL VARIABLES

User CV Name	Stability T/S	Linearity T/S
EFB Bypass	0	0
RFampMode	1	1
Sense Loop	0	0
StabNoGrad	1	0
Lin&Fidlty	0	1
Stab/Grad	4	0
Bandwidth	0	0
BW Cal	0	0

6. Start **[Manual Prescan]**.
7. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".
 Plot Type = **Magnitude**
 Plot Gain = **1**
 Plot Type = **Magnitude**
 Plot Gain = **1**
8. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
9. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
10. Select **[Done]**.
11. Select **[Scan]**. When scan completes, view results by using the Report Manager tool. Record data for each test in Data Sheet 0.7T found in the appendix.


5-5 Coil Tune Preamplifier Check

1. Refer to Illustration 5-9. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
 ILLUSTRATION 5-9

2. Bypass the TNF unit at the rear of the System Cabinet.
3. Refer to Illustration 5-10 for the hardware setup:



COIL TUNE CHECK PREAMPLIFIER LOOP BACK MANUAL PRESCAN SETUP
 ILLUSTRATION 5-10

4. Prepare the system to scan in Body mode per Table 5-9.

TABLE 5-9
BODY SCAN PROTOCOL: BODY COIL TUNE PREAMPLIFIER CHECK TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen] (see Table 5-2 for CV Values):

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **7**<ENTER>.

CV Name: **Sense Loop**, CV Value: **1**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **0**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **0**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **0**<ENTER>.

CV Name: **Bandwidth**, CV Value: **1**<ENTER>.

CV Name: **BW Cal**, CV Value: **1**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>


[Done].

F. **[Research Operations]** **[Download]**.

5. Start [**Manual Prescan**].
6. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".
Plot Type = **Magnitude**
Plot Gain = **1**
Plot Type = **Magnitude**
Plot Gain = **1**
7. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 over range. Do NOT adjust R1 or R2.
8. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
9. Select [**Done**].
10. Select [**Scan**]. This scan data is used by RFT Analysis for frequency compensation with hardware configurations (Body Coil Tune Sense Loop Check) and (Head Coil Tune Sense Loop Check). (Viewing RFT results for this scan not necessary.)


5-6 Preamp Bandwidth Check

1. Refer to Illustration 5-11. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the DIS (disable faults) position.
 - DD switch to the DIS (disable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES DISABLED
 ILLUSTRATION 5-11

2. Bypass the TNF unit at the rear of the System Cabinet.
3. Refer to Illustration 5-12 for hardware setup.



PREAMPLIFIER BANDWIDTH CHECK TPS LOOPBACK MANUAL PRESCAN
 ILLUSTRATION 5-12

4. Prepare the system to scan in Body mode per Table 5-11.

TABLE 5-11

BODY SCAN PROTOCOL: BODY COIL PREAMPLIFIER BANDWIDTH CHECK TEST

Note: This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

A. **[New Pt]**

Id: **geservice**<ENTER>

Name: **rft**

Weight (Lb.): **300**<ENTER>

Set Patient Protocols to **Service**.

At front enclosure, landmark in head area, press **LANDMARK**, then **MOVE TO SCAN**.

B. In the Patient Position Protocol field:

type **o.16.1**<ENTER>(o=Other, 16.1 =series) to load the body protocol

OR select **other** and select protocol **16** (RF Test) and select series **1** (RFT Body).

C. In the Additional Parameters field:

[User CVs Screen] (see Table 5-2 for CV Values):

CV Name: **EFB Bypass**, CV Value: **0**<ENTER>.

CV Name: **RFampMode**, CV Value: **7**<ENTER>.

CV Name: **Sense Loop**, CV Value: **0**<ENTER>.

CV Name: **StabNoGrad**, CV Value: **0**<ENTER>.

CV Name: **Lin&Fidlty**, CV Value: **0**<ENTER>.

CV Name: **Stab/Grad**, CV Value: **0**<ENTER>.

CV Name: **Bandwidth**, CV Value: **1**<ENTER>.

CV Name: **BW Cal**, CV Value: **0**<ENTER>.

[Accept].

D. **[Save Series]**. If necessary --- **[Prepare to Scan]**.

E. **[Research Operations]** **[Setup Params]**.

R1 = 13

R2 = 14

TG = 50

Number of Frames: **4**<ENTER>

Window 1: Frame: **1**<ENTER>

+/- = +

Window 1: Frame: **0** <ENTER>

Window 2: Frame: **2** <ENTER>

+/- = +

Window 2: Frame: **0** <ENTER>

[Done].

F. **[Research Operations]** **[Download]**.

5. Start [**Manual Prescan**].
6. From the "Windows" Menu bar on the Manual Prescan window, select "**Two Windows**".
Plot Type = **Magnitude**
Plot Gain = **1**
Plot Type = **Magnitude**
Plot Gain = **1**
7. Advance TG from 50 to 200, while adjusting variable attenuators A and B to prevent R1 and R2 overrange. Do NOT adjust R1 or R2.
8. With TG at 200, make final adjustments to variable attenuators A and B for R2 signal between 75 & 85%. Do NOT adjust R1 or R2.
9. Select [**Done**].
10. Select [**Scan**]. When scan completes, view results by using the Report Manager tool. Isolate/troubleshoot the bandwidth problem based on the results of this test.


6- RESTORATION CHECKLIST

1. Remove Sense Coils and foam from inside of Body Coil or Head Coil.
2. At the System Cabinet:
 - A. Remove Switchable Attenuator and associated test hardware.
 - B. Enable the TNF (Transient Noise Filter) unit in the rear of the System Cabinet.
 - C. Re-connect all System Cabinet cables to original configuration.
3. At the RF Cabinet:
 - A. Remove associated test hardware.
 - B. Re-connect RF Cabinet cables to original configuration.
4. At the rear pedestal:
 - A. Remove associated test hardware.
 - B. Re-connect Body Preamplifier to original configuration.



Equipment damage possibility. Application of RF energy to TR Switches with Jumpers in TUNE position can cause permanent damage or destruction of components. The Jumpers must be placed back to the scan (normal) position. This JP1 and JP2 position differs for 1.5T and 1.0T Head Coils. Refer to the Head Coil Tuning Sense Loop Test procedure for specific details.

5. **For 1.5T Head Coil Tune Only:** Reinstall rigid coax's and install JP1 and JP2 jumpers on Isolation Boards in the B (scan) position.
6. **For 1.0T Head Coil Tune Only:** Install JP1 and JP2 jumpers in the A (scan) position.
7. **For Body and Head Coil Tune Only:** Remove N shorting caps from I and Q lines and reconnect I and Q lines to the Body Hybrid Splitter.
8. Refer to Illustration 6-1. On the front of the SSM place the DRV FLT switches as follows:
 - TR switch to the **EN** (enable faults) position.
 - DD switch to the **EN** (enable faults) position.



SSM FRONT PANEL TR AND DD FAULT SWITCHES ENABLED
ILLUSTRATION 6-1

7. When finished with RFT, delete all unneeded RFT studies.
8. Perform at least one head or body scan to verify proper system operation.

7- RFT RESULTS

7-1 RFT Analysis

RFT data analysis begins immediately after the scan is completed. The following is displayed on the display plasma while the data analysis is occurring:

Note

Only those tests which were selected during scan prescription will be processed.

Running RF Transmit Test Analysis

```
*****  
*                               RFT Analysis Tool                               *  
*****
```

```
RF Stability:Slice 1, echo 1  
RF Stability:Slice 2, echo 1  
RF Stability:Slice 3, echo 1  
RF Stability:Slice 4, echo 1  
RF Linearity:Slice 1, echo 2  
RF Linearity:Slice 2, echo 2  
RF Droop:Slice 2, echo 1  
RF Fidelity:Slice 2, echo 3  
RF Fidelity:Slice 2, echo 4  
RF Bandwidth:Slice 5, echo 1
```

Note: RF Bandwidth appears only when Bandwidth = 1.

```
***** RFT Analysis Completed *****
```

Press any key to quit-->Enter any key to quit.

7-2 Report Program

RFT results (both text and graphics) can be viewed using the Report program from the Signa Tools menu. Refer to Table 7-1. This section describes how to initiate the Report program. Remaining sections provide examples of all reported parameters for RFT results. Record data for each test in Data Sheet 1.5T or Data Sheet 1.0T found in Appendix A.

Note

In Release 8.2 (and later), the RFT data sheet can automatically be generated. Refer to the Report Manager procedure.

TABLE 7-1
 RESULT FILE ORGANIZATION

Group	Slice	Echo	Title	Plots
1	1	1	Low Power Magnitude Stability - - No Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
2	1	1	Low Power Phase Stability - - No Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
3	2	1	High Power Magnitude Stability - - No Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
4	2	1	High Power Phase Stability - - No Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
5	3	1	Low Power Magnitude Stability - - With Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
6	3	1	Low Power Phase Stability - - With Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation

TABLE 7-1 - CONTINUED
 RESULT FILE ORGANIZATION

Group	Slice	Echo	Title	Plots
7	4	1	High Power Magnitude Stability - - With Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
8	4	1	High Power Phase Stability - - With Gradients	1 Sample deviation all (max) 2 Sample deviation all (min) 3 View means 4 Min sample deviation 5 Max sample deviation
9	1 & 2	2	Linearity -61 dB to 0 dB	1 Magnitude Error 2 Phase Error
10	2	1	Magnitude Droop	1 Raw Droop 2 Early Droop 3 Average Droop 4 Late Droop
11	2	3	Sinc 2 Pulse Fidelity	1 Half Scale Odd Views 2 Full Scale Odd Views 3 Half Scale Even Views 4 Full Scale Even Views 5 Ideal Slice
12	2	4	Sinc 8 Pulse Fidelity	1 Half Scale Odd Views 2 Full Scale Odd Views 3 Half Scale Even Views 4 Full Scale Even Views 5 Ideal Slice
13	5	1	Bandwidth	1 Raw Magnitude vs. Frequency 2 Corrected Magnitude vs. Frequency
14	5	1	Coil Tune Check	1 Raw Magnitude vs Frequency 2 Corrected Magnitude vs. Frequency

7-3 RFT Header Data

Refer to Table 7-2 for definition of legal value entries for all RFT Report header parameters:

TABLE 7-2
RFT HEADER PARAMETERS

Parameter	Legal Values	Source
Sitename	Site Name	Raw Header
Usn	Unique System Number (GE Cares Issued)	Config File
Srvconfig	Date/Time Srv Config File Last Changed	Srv Config File ("sysconfig")
Receiver	# Of Tps Rcvr Bds (4-4); # Channels/Bd (1) Starting Rec., Ending Rec, Port Enable	Srv Config File
Xmtrfcoil	Body, Head	Raw Header
Rcvrfcoil	Body, Head, <Surface Coil Name>	Raw Header
Freq	Magnet Frequency	Raw Header
Time	YY/MM/DD HH:MM:SS	Raw Header
Baserun	Base Run Number Of Scan	Raw Header
Configcode	_xxxxxxxxxx	Raw Header
Softrev	Software Revision	"mrswrev" Script
Rcvcoilgain	Calibration Value (Typ. 4-10)	Config File
R1	4-7	Raw Header
R2	4-30	Raw Header
TG	0-200	Raw Header
C1	Exam Description	Raw Header
C2	Series Description	Raw Header

Refer to the following for a sample header:

```
17JE1325.RFT/0 Header Info RDF/GRP Revision: 1.0 /41
=====
SITENAME      = BAY 12
USN           = B15B0C0
MLN           = 9999
SRVCONFIG     = YY/MM/DD HH/MM/SS
EXCITER       = 000
RECEIVER      = S0/E0/PE0
XMTRFCOIL     = BODY
RCVRF0IL     = BODY
FREQ          = 42663520 Hz
TIME          = 91/7/27 07:40:00
BASERUN       = 44032
CONFIGCODE    = _0201240000
SOFTREV       = 5.3
NUCLIDE       = 000
HEADERCODE    = 0x00000000
RCVCOILGAIN  = 6.090  R1 = 13  R2 = 14  TG = 200
----- Exam description -----
NO COMMENT
----- Series description -----
Body,Ax,2D,Spin Echo,None
-----
```

7-4 RFT Sample Test Outputs

Refer to the following sections for sample RFT data for each test mode.


7-4-1 Low Power Magnitude (No Gradients) Stability Data

```
xxxxxxx.RFT/1 LO PWR MAGNITUDE STABILITY - NO GRADIENTS
=====
All values below are percent of peak pulse amplitude.
-----
max-min(pp)   + stdev   - stdev   + area   - area
    0.20      0.01     0.01     14.46   -14.99
-----
max sample column 157      p-p      stdev
min sample column  32      0.16     0.02
view means          0.07     0.01
-----
```

7-4-2 Low Power Phase (No Gradients) Stability Data

```

xxxxxxxx.RFT/2  LO PWR PHASE STABILITY - NO GRADIENTS
=====
All values below are phase deviation in degrees.
-----
max-min(pp)    + stdev      - stdev      + area      - area
   0.47         0.02         0.03         18.79       -28.01
-----
max sample column  5          p-p          stdev
min sample column  5          0.47         0.09
view means         0.12         0.02
-----
  
```



RFT LOW POWER PHASE STABILITY (NO GRADIENTS) SAMPLE PLOTS
 ILLUSTRATION 7-1

7-4-3 High Power Magnitude (No Gradients) Stability Data

xxxxxxxx.RFT/3 HI PWR MAGNITUDE STABILITY - NO GRADIENTS


=====

All values below are percent of peak pulse amplitude.

max-min(pp)	+ stdev	- stdev	+ area	- area
0.28	0.01	0.03	10.26	-37.37

	p-p	stdev
max sample column 10	0.16	0.02
min sample column 180	0.21	0.04
view means	0.17	0.03

Plot below may be incorrect-refer to Illustration 7-1.




**HIGH POWER MAGNITUDE STABILITY (NO GRADIENTS) SAMPLE PLOTS
 ILLUSTRATION 7-2**

7-4-4 High Power Phase (No Gradients) Stability Data

```

xxxxxxxx.RFT/4  HI PWR PHASE STABILITY - NO GRADIENTS
=====
All values below are phase deviation in degrees.
-----
max-min(pp)    + stdev      - stdev      + area      - area
    0.50        0.04        0.03        49.85      -32.21
-----
max sample column 5          p-p          stdev
min sample column 206       0.37        0.08
view means          0.40        0.05
                    0.29        0.04
-----
    
```

Plot below may be incorrect-refer to Illustration 7-1.



HIGH POWER PHASE STABILITY (NO GRADIENTS) SAMPLE PLOTS
 ILLUSTRATION 7-3

7-4-5 Low Power Magnitude (With Gradients) Stability Data


xxxxxxxx.RFT/5 LO PWR MAGNITUDE STABILITY WITH GRADIENTS

=====
 All values below are percent of peak pulse amplitude.

max-min(pp)	+ stdev	- stdev	+ area	- area
0.41	0.02	0.06	13.21	-32.16

	p-p	stdev
max sample column 212	0.15	0.03
min sample column 33	0.35	0.05
view means	0.10	0.02

Plot below may be incorrect-refer to Illustration 7-1.




LOW POWER MAGNITUDE STABILITY (WITH GRADIENTS) SAMPLE PLOTS
 ILLUSTRATION 7-4

7-4-6 Low Power Phase (With Gradients) Stability Data

```

xxxxxxxx.RFT/6  LO PWR PHASE STABILITY WITH GRADIENTS
=====
All values below are phase deviation in degrees.
-----
max-min(pp)    + stdev      - stdev      + area      - area
   0.42         0.02         0.03         15.76       -33.28
-----
max sample column 6          p-p          stdev
min sample column 250       0.41         0.06
view means          0.34         0.06
                    0.12         0.02
-----
  
```

Plot below may be incorrect-refer to Illustration 7-1.




LOW POWER PHASE STABILITY (WITH GRADIENTS) SAMPLE PLOTS
 ILLUSTRATION 7-5

7-4-8 High Power Phase (With Gradients) Stability Data

```

xxxxxxx.RFT/8  HI PWR PHASE STABILITY WITH GRADIENTS
=====
All values below are phase deviation in degrees.
-----
max-min(pp)    + stdev      - stdev      + area      - area
    0.45        0.05        0.01        57.23      -16.91
-----
max sample column 7          p-p          stdev
min sample column 219       0.42        0.08
view means           0.34        0.06
view means           0.25        0.06
-----
    
```

Plot below may be incorrect-refer to Illustration 7-1.



HIGH POWER PHASE STABILITY (WITH GRADIENTS) SAMPLE PLOTS
 ILLUSTRATION 7-7

7-4-9 Linearity Data

xxxxxxxx.RFT/9 LINEARITY -61 to 0 dB
 =====
 Carrier Leakage (dB below full scale modulation)

N/A


 Magnitude Linearity Error (dB)

max	Min	P-P	Std Dev
0.23	-0.02	0.25	0.06

 Phase Linearity Error (degrees)

max	Min	P-P	Std Dev
2.41	-3.75	6.17	2.09

Plot below may be incorrect for 8X.



LINEARITY SAMPLE PLOT
 ILLUSTRATION 7-8


7-4-10 Magnitude Droop Data

```

xxxxxxxx.RFT/10  MAGNITUDE DROOP
=====
Peak Amplitude  27260.2 counts
-----
Droop Rates in dB/millisecond

Early          0.0820
Average        0.0030
Late           -0.0000
-----
  
```

Plot below may be incorrect for 8X.



MAGNITUDE DROOP SAMPLE PLOT
 ILLUSTRATION 7-9

L1775A

7-4-11 Sinc 2 Pulse Fidelity Data

```

xxxxxxxx.RFT/11      SINC 2 PULSE FIDELITY
=====
Positive Modulation


Full Scale
  Slice Amplitude      121.71 % of ideal slice
  Total Error          303.27 dB x samples > -45 dB
  Pulse Amplitude     81.15 % of ideal pulse

Half Scale
  Slice Amplitude      58.02 % of ideal slice
  Total Error          290.59 dB x samples > -45 dB
-----
Negative Modulation

Full Scale
  Slice Amplitude      121.04 % of ideal slice
  Total Error          318.04 dB x samples > -45 dB
  Pulse Amplitude     80.98 % of ideal pulse

Half Scale
  Slice Amplitude      57.23 % of ideal slice
  Total Error          347.59 dB x samples > -45 dB
-----
  
```

Plot below may be incorrect for 8X.



SINC 2 PULSE FIDELITY SAMPLE PLOT
 ILLUSTRATION 7-10

L1776A

7-4-12 Sinc 8 Pulse Fidelity Data

xxxxxxxx.RFT/12 SINC 8 PULSE FIDELITY

=====
 Positive Modulation

Full Scale

Slice Amplitude 30.25 % of ideal slice
 Total Error 1655.74 dB x samples > -45 dB
 Pulse Amplitude 80.45 % of ideal pulse

Half Scale

Slice Amplitude 14.41 % of ideal slice
 Total Error 2030.24 dB x samples > -45 dB

 Negative Modulation


Full Scale

Slice Amplitude 30.19 % of ideal slice
 Total Error 1683.67 dB x samples > -45 dB
 Pulse Amplitude 80.40 % of ideal pulse

Half Scale

Slice Amplitude 14.35 % of ideal slice
 Total Error 2127.08 dB x samples > -45 dB

Plot below may be incorrect for 8X.




SINC 8 PULSE FIDELITY SAMPLE PLOT
 ILLUSTRATION 7-11

L1777A

7-4-13 Bandwidth Data

```
xxxxxxxxx.RFT/13      BANDWIDTH
=====
Center Frequency      -0.25 dB Bandwidth
    63.777 MHz                0.346 MHz
-----
Low -0.25 dB freq    High -0.25 dB freq
    63.612 MHz                63.958 MHz
-----
Sweep start freq     Sweep end freq
    63.360 MHz                64.360 MHz
-----
Cal File Data Acquired
-----
```

Plot below may be incorrect for 8X.



RFT BANDWIDTH SAMPLE PLOT
ILLUSTRATION 7-12


L1778A

7-4-14 Coil Tune Check Data

```

xxxxxxxxx.RFT/14    COIL TUNE CHECK
=====
Center Frequency          -3.0 dB Bandwidth
    63.777 MHz                0.346 MHz
-----
Low -3.0 dB freq         High -3.0 dB freq
    63.612 MHz                63.958 MHz
-----
Sweep start freq         Sweep end freq
    63.360 MHz                64.360 MHz
-----
Cal File Data Acquired
-----
    
```

Plot below may be incorrect for 8X.



COIL TUNE CHECK SAMPLE PLOT
 ILLUSTRATION 7-13

8L1779A

7-4-15 Config File Data

```
xxxxxxxx.RFT/15/15    SERVICE CONFIG FILE
=====
magnetSerNum          = "R244"
magnetType             = "GELCC-15"
fieldStrength         = "1.5T"
room                  shielding (Yes/No) = "No"
magniShield           = (Yes/No) = "No"
shieldCooler          = type (0,1,2,3,4,5,6) = "4"
resistive              shim power supply (0,1,2,3,4,5,6) = "0"
RFAmpType = (ErbtecTube, ErbtecSolidState, TELI Solid State) = "4"
spectroAmp = (none, ENI) = "None"
line CondType         (0,1,2) = "2"
lineFreq              (50, 60) = "60"
SRMode                (120,77,20,17) = "120"
ThisLineNotUsed      = ""
ThisLineNotUsed      = ""
ThisLineNotUsed      = ""
-----
```

8- RFT THEORY

8-1 RFT Overview

RFT is an RF loop-back test that will quickly evaluate stability, magnitude and phase linearity, droop, pulse fidelity, and bandwidth of the RF transmit chain in a Signa system. This is accomplished by playing out four different RF pulses and collecting each one as a separate "echo". The first pulse has a wide flat top and is used for stability, droop and bandwidth testing. The second pulse has a narrower flat top and its amplitude increments during the scan for magnitude and phase linearity testing. The third pulse is a sinc 2 and the fourth pulse is a sinc 8. The amplitude of the third and fourth pulses is incremented like the second pulse for pulse fidelity analysis.

All tests are run at the following frequencies: (Actual magnet center frequency will not be used.)

1.5T - 63.86 MHz

1.0T - 42.57 MHz

For the bandwidth test, provision is made to characterize the Exciter/Receiver combination so that when the RF amp is tested the results can be corrected for errors in the Exciter/Receiver. When testing bandwidth using a sense loop in one of the scanner's transmit coils, the results will likewise be corrected for errors in the Exciter/RF Amplifier combination. User CVs are provided so the analysis tool can be told the exact status of the loop back configuration.

RFT Hardware

- A switchable attenuator controlled by the PSD through the DAB is required to extend the dynamic range of the linearity and stability tests. Linearity and stability data are collected twice, first with the transmit path attenuator set to 30 dB (IN) and again with it set to 0 dB (OUT). The receive path attenuator is 30 dB when the transmit path attenuator is 0 dB and vice versa. Stability data is collected twice more with a gradient pulse applied during the RF pulse. The switchable attenuator is used as before.


Testing Sequence

- This document is structured for four basic modes of running RFT: 1) baseline (reference) data collection, 2) coil tune check, 3) stability troubleshooting, and 4) pulse fidelity troubleshooting. Because the number of possible hardware and software configurations is large, flowcharts are provided for each mode. The flowcharts direct you to the appropriate tables and diagrams for software and hardware setup, prescan, and scan references.

8-2 Data Collection

First slice

The 128 views are collected with the switchable attenuator IN (i.e., transmit path set to 30 dB). Echo 2, 3, and 4 pulse amplitudes are incremented as specified in Illustration 8-1.




SLICE 1 & 2 DATA COLLECTION
 ILLUSTRATION 8-1

Second slice

The 128 views are collected with the switchable attenuator OUT (i.e., transmit path set to 0 dB). Echo 2, 3, and 4 pulse amplitudes are incremented as specified in Illustration 8-1.

Third slice

The 128 views are collected with the switchable attenuator IN (i.e., transmit path set to 30 dB). Echo 2, 3, and 4 pulse amplitudes are set to zero with modulator unblank, EFB enable and Data In Windows disabled as specified in Illustration 8-2.



SLICE 3, 4, & 5 DATA COLLECTION
 ILLUSTRATION 8-2

Fourth slice

The 128 views are collected with the switchable attenuator OUT (i.e., transmit path set to 0 dB). Echo 2, 3, and 4 pulse amplitudes are set to zero with modulator unblank, EFB enable and Data In Windows disabled as specified in Illustration 8-2.

Fifth slice

The switchable attenuator is OUT (i.e. transmit path set to 0 dB) and operating frequency is set to 500 kHz below 63.86 MHz (42.57 or 21.89 MHz). Transmit and receive frequency are incremented (1000/127) kHz per view while the echo 1 pulse is played out at full scale. The pulse amplitude for echoes 2, 3, and 4 is zero and the exciter and RF amplifier are not unblanked during the 128 views of the third slice. See Illustration 8-2.

8-3 Magnitude Stability Analysis (groups 1, 3, 5 and 7)

A magnitude matrix is calculated from echo 1 raw data. Each element in the normalization vector (**nvc**) contains the mean of the corresponding column in the magnitude matrix. Peak magnitude (**pkmag**) is the mean of elements 31 through 37 in the normalization vector. Each element in the sample max vector (**smaxc**) contains the maximum of the corresponding column in the normalized magnitude matrix. (**cmax**) is element number (**c**) for which smaxc is maximum. Each element in the sample min vector (**sminc**) contains the minimum of the corresponding column in the normalized magnitude matrix. (**cmin**) is element number (**c**) for which sminc is minimum. Each element in the view means vector (**vmeanr**) contains the mean of the corresponding row in the normalized magnitude matrix.

The test provides plots of "Sample Deviation All" (max% & min%), "View Means" (%), "Max Sample Deviation" (%), and "Min Sample Deviation" (%). Report data file summary:

- Max-Min (p-p)
- +Stdev
- -Stdev
- +Area
- -Area
- Max Sample Column
- Min Sample Column
- View Means

8-4 Phase Stability Analysis (groups 2, 4, 6 & 8)


An unwrapped phase matrix is calculated from echo 1 raw data. Each element in the normalization vector (**nvc**) contains the mean of the corresponding column in the unwrapped phase matrix. Each element in the sample max vector (**smaxc**) contains the maximum of the corresponding column in the normalized phase matrix. (**cmax**) is element number (**c**) for which smaxc is maximum. Each element in the sample min vector (**sminc**) contains the minimum of the corresponding column in the normalized phase matrix. (**cmin**) is element number (**c**) for which sminc is minimum. Each element in the view means vector (**vmeanr**) contains the mean of the corresponding row in the normalized phase matrix.

The test provides plots of "Sample Deviation All" (max% & min%), "View Means" (%), "Max Sample Deviation" (%), and "Min Sample Deviation" (%). Report data file summary:


- Max-Min (p-p)
- +Stdev
- -Stdev
- +Area
- -Area
- Max Sample Column
- Min Sample Column
- View Means

8-5 Linearity And Carrier Leakage Analysis (group 9)

Data for Linearity and Carrier Leakage Analysis is collected in passes (slices) 1 and 2 of the scan. The second echo is analyzed for linearity. Polarity of the RF modulation is alternated on every other view and 4 views are taken at each magnitude. Therefore, there are 2 positive pulses and 2 negative pulses for each amplitude. Amplitude is incremented in 1.0 dB steps beginning at 31.0 dB below full scale and ending at full scale for a total of 32 different amplitudes in each pass. The switchable attenuator inserts 30 dB between output of the exciter and input of the RF amplifier or in the loopback path between output of the RF amplifier and input of the receiver. See Illustration 8-3 and Illustration 8-4. Since pulse amplitude is incremented through 32 dB in each pass and the attenuation change is 30 dB between passes, there is 2 dB of overlap in RF amplifier output between passes 1 and 2. This overlap is used to smoothly splice together linearity data for the two passes yielding a total dynamic range for the linearity test of 62 dB.



PASS 1 & 3 SWITCHABLE ATTENUATOR CONDITION
ILLUSTRATION 8-3



PASS 2, 4, & 5 SWITCHABLE ATTENUATOR CONDITION
ILLUSTRATION 8-4

The test provides plots of "Magnitude Error" (dB), and "Phase Error" (deg). Report data file summary:

- Carrier Leakage
- Magnitude Linearity Error (Max)
- Magnitude Linearity Error (Min)
- Magnitude Linearity Error (p-p)
- Magnitude Linearity Error (Std Dev)
- Phase Linearity (Max)
- Phase Linearity (Min)
- Phase Linearity (p-p)
- Phase Linearity (Std Dev)

8-6 Magnitude Droop Analysis (group 10)

Data for this test is collected in pass (slice) 2 of the scan. The first echo is analyzed for magnitude droop.

The test provides plots of "Raw Droop", "Early Droop", "Average Droop" and "Late Droop".

Report data file summary:

- Peak Amplitude
- Early Droop Rate
- Average Droop Rate
- Late Droop Rate

8-7 Sinc 2 Pulse Fidelity Analysis (group 11) (TPS Loopback only)

Data for this test is collected in pass (slice) 2 of the scan. The third echo is analyzed for pulse fidelity.

The test provides plots of "+Half Scale Slice Error" (dB), "+Full Scale Slice Error" (%), "-Half Scale Slice Error" (%) and "-Full Scale Slice Error" (dB). Report data file summary:

- + Full Scale Slice Amplitude
- + Full Scale Total Error
- + Full Scale Pulse Amplitude
- + Half Scale Slice Amplitude
- + Half Scale Total Error
- - Full Scale Slice Amplitude
- - Full Scale Total Error
- - Full Scale Pulse Amplitude
- - Half Scale Slice Amplitude
- - Half Scale Total Error

8-8 Sinc 8 Pulse Fidelity Analysis (group 12) (TPS Loopback only)

Data for this test is collected in pass (slice) 2 of the scan. The fourth echo is analyzed for pulse fidelity.

The test provides plots of "+Half Scale Slice Error" (dB), "+Full Scale Slice Error" (%), "-Half Scale Slice Error" (%) and "-Full Scale Slice Error" (dB). The Report data file summary is:

- + Full Scale Slice Amplitude
- + Full Scale Total Error
- + Full Scale Pulse Amplitude
- + Half Scale Slice Amplitude
- + Half Scale Total Error
- - Full Scale Slice Amplitude
- - Full Scale Total Error
- - Full Scale Pulse Amplitude
- - Half Scale Slice Amplitude
- - Half Scale Total Error


8-9 Magnitude Bandwidth Analysis (group 13 & 14)

Data for this test is collected in pass (slice) 5 of the scan. The first echo is analyzed for magnitude bandwidth.

Reference files are stored in /usr/g/service/data. They are required for the following conditions:

TPS Reference File (Group 13):

- The bandwidth test is run with output of the TPS looped back to the input of the receiver via the switchable attenuator as shown in Illustration 8-5. See Table 8-1 to determine when a copy of the filtered plot data will be written to the file /usr/g/service/data/tpsbw.cal.



L1725A

TPS LOOPBACK SETUP
ILLUSTRATION 8-5

RF Amp Body Mode Reference File (Group 13):

- The bandwidth test will be run with output of the RF amplifier looped back through suitable attenuation to the input of the receiver via the switchable attenuator as shown in (Body Dummy Load, EFB Bypassed). See Table 8-1 to determine when a copy of the filtered plot data (prior to TPS bandwidth correction) will be written to the file /usr/g/service/data/bodybw.cal.

TABLE 8-1
 BW CAL FILE LOGIC

OPUSER 1	OPUSER 2	OPUSER 7	READ CAL FILE	WRITE CAL FILE	RESULT FILE NUMBER	DBDOW N
0	0	0	none	none	13	-.5
0	0	1	none	tpsbw.cal	13	-.5
1	0	0	tpsbw.cal	none	13	-.5
2	0	0	tpsbw.cal	none	13	-.5
2	0	1	tpsbw.cal	headbw.cal	13	-.5
2	1	0	headbw.cal	none	13	-.5
3	0	0	tpsbw.cal	none	13	-.5
4	0	0	tpsbw.cal	none	13	-.5
4	0	1	tpsbw.cal	bodybw.cal	13	-.5
4	1	0	bodybw.cal	none	13	-.5
7	0	0	tpsbw.cal	none	13	-.5
7	0	1	tpsbw.cal	preampbw.cal	13	-.5
7	1	0	preampbw.cal	none	14	-.3

RF Amp Head Mode Reference File (Group 13):

- The bandwidth test will be run with output of the RF amplifier looped back through suitable attenuation to the input of the receiver via the switchable attenuator as shown in (Head Dummy Load, EFB Bypassed). See Table 8-1 to determine when a copy of the filtered plot data (prior to TPS bandwidth correction) will be written to the file /usr/g/service/data/headbw.cal.

Coil Tune Check Preamp Mode Reference File (Group 14):

- The bandwidth test will be run with output of the body preamplifier looped back through suitable attenuation to the input of the receiver as shown in Body tests (Coil Tune Body Preamplifier) & (Body Coil Tune Sense Loop), and Head tests (Coil Tune Head Preamplifier) & (Head Coil Tune Sense Loop). See Table 8-1 to determine when a copy of the filtered plot data (prior to TPS bandwidth correction) will be written to the file /usr/g/service/data/preampbw.cal.
- The test provides plots of "Raw" (MHz), and "Corrected" (MHz). Report data file summary:
 - Center Frequency
 - -0.5dB Bandwidth
 - Low -0.5dB Frequency
 - High -0.5dB Frequency
 - Sweep Start Frequency
 - Sweep End Frequency
 - Cal File Data Acquired

APPENDIX A - DATA SHEET

0.7T RFT BODY/HEAD DATA

Site:		Name:		Date:
RFT File Number:				
Coil:	<input type="checkbox"/> Body <input type="checkbox"/> Head	<input type="checkbox"/> Body <input type="checkbox"/> Head	Body/Head Spec.	
LO PWR MAGNITUDE STABILITY - NO GRADIENTS			5.0%*(1.0% no EFB) (MAX - MIN (PP) %)	
LO PWR PHASE STABILITY - NO GRADIENTS			10.0°* (1.5° no EFB) (MAX - MIN (PP) °)	
HI PWR MAGNITUDE STABILITY - NO GRADIENTS			1.0% (MAX - MIN (PP) %)	
HI PWR PHASE STABILITY - NO GRADIENTS			1.5° (MAX - MIN (PP) %)	
LO PWR MAGNITUDE STABILITY - WITH GRADIENTS			1.0% (MAX - MIN (PP) %)	
LO PWR PHASE STABILITY - WITH GRADIENTS			1.5° (MAX - MIN (PP) %)	
HI PWR MAGNITUDE STABILITY - WITH GRADIENTS			1.0% (MAX - MIN (PP) %)	
HI PWR PHASE STABILITY - WITH GRADIENTS			1.5° (MAX - MIN (PP) %)	
MAGNITUDE LINEARITY -61 to 0 dB			2.5 dB p-p (MAX (dB))	
PHASE LINEARITY -61 to 0 dB			12 %* (MAX - MIN (PP) %)	
MAGNITUDE DROOP			Not Available	
BANDWIDTH			.5 dB down @ +/- 100 kHz from CF (MHz)	
COIL TUNE CHECK			BODY: CF (MHz) +/- 50 kHz HEAD: NO SPEC	

Note that the numbers for Head Dummy Load (less RF drive) will have phase linearity numbers about 40% lower than those in the Body mode.

APPENDIX B - CV DESCRIPTIONS

Name	Meaning	Controls	Allowable Values	Conditions required for values other than zero
EFB Bypass (opuser0)	0= No 1= Yes	Sets allowable values for RFampMode (opuser1) & Bandwidth (opuser6)	[0.0 to 1.0]	None NOTE: Set to 0 for RF/PDU Cabinet.
RFampMode (opuser1)	0= Exciter 1= Head Dummy Load 2 = Head coil load 3= Body Dummy Load 4= Body coil load 5= Test 6= SSdriver 7= Coil tune check	RF Amp Operating Mode 0= Standby 1= Head Mode 2 = Head mode 3= Body Mode 4= Body mode 5= Test mode 6= Head Mode 7= Standby	[0.0 to 7.0]	If RFampMode = 5 or 6 then EFB bypass must = 1. If RFampMode = 1, 2 or 6, then op coil must = 1 (Head coil) If RFampMode = 3 or 4, then op coil must = 2 (Body coil) Note: opcoil is coil selection
Sense Loop (opuser2)	0= No 1= Yes	Sets allowable valuefor BW Cal (opuser7)	[0.0 to 1.0]	RFampMode = 2, 4 or 7 only
StabNoGrad (opuser3)	0 = Tests Disabled 1= Yes	0 = No Stability Analysis 1 = Do Stability Analysis Note: Slices 1 & 2, echo 1 will always be played out even if StabNoGrad = 0.	[0.0 to 1.0]	None
Lin&Fidlty (opuser4)	0 = Test Disabled 1 = Lin only 2 = Lin & Fidelity	0 = Slices 1 & 2, echo 1 only 1 = Slices 1 & 2, echoes 1 & 2 2 = Slices 1 & 2, 4 echoes	[0.0 to 2.0]	If Lin&Fidlty = 2 then RFampMode must = 0 (Exciter)
Stab/Grad (opuser5)	0 = Test Disabled 1 = X Gradient Only 2 = Y Gradient Only 3 = Z Gradient Only	0 = Slices 3 & 4 not done 1 = do slices 3 & 4 w/ X grad. 2 = do slices 3 & 4 w/ Y grad 3 = do slices 3 & 4 w/ Z grad	[0.0 to 4.0]	None

	4= All Gradients	4 = do slices 3 & 4 w/ all grad. Note: Slices 3 & 4 will be played out but not analyzed if Stab/Grad = 0 and Bancwidth = 1.		
Bandwidth (opuser6)	0 = Test Disabled 1 = Test Enabled	0 = slice 5 not done 1 = Do slice 5 Sets allowable value for BW Cal	[0.0 to 1.0]	RFampMode = 0 or 7 OR RFampMode = 1, 2, 3 or 4 AND EFB Bypass = 1
BW Cal (opuser7)	0 = Disabled 1 = Enabled	Analysis Use Only	[0.0 to 1.0]	Bandwidth = 1 AND RFampMode = 0, 2, 4 or 7 AND Sense Loop = 0

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
A	Dec 1, 1999	Resa Lambert	Preliminary version.
0			Initial version.